

# AssureCare LLC

# 170.315(b) (10) Electronic Health Information

# Export

Version 1.0.0.1

Last updated: November 29, 2022

The material presented in this document is intended for use of iPatientCare clients only and may not be reproduced in any form, by any method, for any purpose without the expressed permission of AssureCare, LLC.

# About Guide

## Revision History

Version	Date	Comment
1.0.0.0	September 15, 2022	Original Document
1.0.0.1	November 29, 2022	Included the Product Name and Version section.

## Conventions

Before using this document, it is important to understand the typographical conventions used to identify and describe information.

Description	Appearance
Field Names, Keys you press and Buttons/Icons you click is shown in Open Sans Semibold.	e.g., Press <b>Enter</b> . e.g., Click <b>Ok</b> to continue. e.g., Select <b>Provider</b> from the drop-down list.
Dialog Box and Application Window Titles are shown in Open Sans italics.	e.g., The <i>Preview</i> popup appears.
Cross-References to section and hyperlink are shown in Open Sans Semibold.	e.g., Refer to <b><a href="#">Conventions Used</a></b> section.
Warning and Alert messages are shown in Italic Quoted Text	e.g., " <i>Warning and Alert messages</i> "
Assumption and example are in Italic	e.g., <i>For example</i>

## Support

You can call, fax or e-mail iPatientCare Support:

Phone: 732.607.2400

E-mail: [support@ipatientcare.com](mailto:support@ipatientcare.com)

For support, you may also generate a ticket from iPatientCare Support Portal:

<http://support.ipatientcare.com>

## Product Name and Version:

iPatientCare 18.0

iPatientCare 22.5

## 170.315(b) (10) Electronic Health Information Export

The applications provide export of electronic health information (EHI) for a single patient as well as for patient population in the following formats:

1. **CCDA:** iPatientCare supports bulk export of HL7 CCDAs which comply to [United States Core Data for Interoperability \(USCDI\)](#), Version 1 requirements.

The specifications for the CCDA can be obtained from the HL7 website:

- [Health Level 7 \(HL7\) Implementation Guide for CDA Release 2: IHE Health Story Consolidation, DSTU Release 1.1 \(US Realm\) Draft Standard for Trial Use July 2012](#)
  - [HL7® Implementation Guide for CDA® Release 2: Consolidated CDA Templates for Clinical Notes \(US Realm\), Draft Standard for Trial Use Release 2.1 August 2015, June 2019 \(with Errata\)](#)
  - [HL7® CDA R2 IG: C-CDA Templates for Clinical Notes R2.1 Companion Guide, Release 2, October 2019, IBR approved for § 170.205\(a\)\(5\)](#)
2. **FHIR:** iPatientCare also supports [HL7® Version 4.0.1 FHIR® Release 4, October 30, 2019](#), [FHIR® US Core Implementation Guide STU V3.1.1](#), [HL7 FHIR Bulk Data export](#). Please refer [170.315\(g\)\(10\) SmartOnFHIR API Documentation.PDF](#) for the API Specifications..