

iPatientCare

SMART-on-FHIR

API Authentication and Access Guide – 1.0.0.3

Last updated: November 12, 2025

About Guide

Revision History

Version	Date	Comment
1.0.0.0	September 15, 2022	Original Document
1.0.0.1	November 29, 2022	Updated the Production Endpoint URLs and content in the Access FHIR API and Terms of Use sections.
1.0.0.2	February 27, 2023	Added version 23.01
1.0.0.3	November 12, 2025	Updated for HL7® FHIR® US Core Implementation Guide STU 6.1.0 , 170.213(b) United States Core Data Interoperability (USCDI), October 2022 Errata, Version (v3) and 170.215(c)(2) HL7® SMART App Launch Implementation Guide Release 2.0.0

Conventions

Before using this document, it is important to understand the typographical conventions used to identify and describe information.

Description	Appearance
Field Names, Keys you press and Buttons/Icons you click is shown in Open Sans Semibold.	e.g., Press Enter. e.g., Click Ok to continue. e.g., Select Provider from the drop-down list.
Dialog Box and Application Window Titles are shown in Open Sans italics.	e.g., The <i>Preview</i> popup appears.
Cross-References to section and hyperlink are shown in Open Sans Semibold.	e.g., Refer to Conventions Used section.
Warning and Alert messages are shown in Italic Quoted Text	e.g., <i>"Warning and Alert messages"</i>
Assumption and example are in Italic	e.g., <i>For example</i>

Support

You can call, fax or e-mail iPatientCare Support:

Phone: 732.607.2400

E-mail: support@ipatientcare.com

For support, you may also generate a ticket from iPatientCare Support Portal:

<http://support.ipatientcare.com>

Contents

FHIR Server Installation.....	6
Product Name and Version	6
Introduction.....	6
FHIR Endpoints.....	6
Production Endpoint.....	6
Test Server Endpoint.....	6
Types of Smart-On-FHIR applications.....	6
Types of Authentications supported.....	7
Steps for Smart App Launch	7
Client APP Registration	7
Obtain authorization code.....	8
Obtain access token	12
Symmetric Authentication Request.....	13
Access FHIR API.....	16
Request : Patient.....	17
Request : AllergyIntolerance	25
Request : CarePlan	28
Request : CareTeam.....	33
Request: Problems / Health Concern	37
Request: Implantable Device	42
Request: Diagnostic Report/ Clinical Notes/ DocumentReference.....	45
Request: Laboratory Result Observation.....	49
Request: Goal	53
Request: Immunization	55
Request: Medication	59
Request: Smoking Status.....	60
Request: Procedure.....	65
Request: Provenance.....	67
Request: Vital Signs	74
Request: Organization.....	87

Request: DocumentReference	89
Request: ServiceRequest.....	93
Request: Coverage	96
Request: MedicationDispense	101
Request: Specimen	105
Request: RelatedPerson	107
Request: Observation.....	109
Profile audience and scope	114
Security and Privacy Considerations.....	115
App Protection	115
Terms of Use	116

FHIR Server Installation

Product Name and Version

iPatientCare 18.0

iPatientCare 22.5

iPatientCare 23.0

Introduction

In this documentation we will list all the steps required to access protected health information based on the documentation of an open API. The API makes health information broadly available using FHIR®, a set of clinical interoperability resources under the umbrella of the HL7 standards organization. FHIR is based on common web standards can be reached through a RESTful protocol in which each FHIR resource has a known URL. This document describes OAuth 2.0 for client applications to authorize, authenticate, and integrate with FHIR-based data systems.

We follow Smart-On-FHIR authentication (OAuth2) and authorization (Ref: <https://hl7.org/fhir/smart-app-launch/index.html>)

FHIR Endpoints

Note: Production and Test Server endpoints are subject to change.

Production Endpoint

- FHIR : <https://fhir.iPatientCare.net:9443/fhirserver/fhir/>
- Authentication Server: <https://oauth.iPatientCare.net/>

Test Server Endpoint

- FHIR : <https://uat-fhir.ipatientcare.net:9443/fhirserver/fhir/>
- Authentication Server: <https://uat-oauth.ipatientcare.net/>

Types of Smart-On-FHIR applications

- Standalone App : SMART on FHIR confidential client with a patient context, refresh token, and OpenID Connect (OIDC) identity token.
- EHR Embedded app: Demonstrate the ability to perform an EHR launch to a SMART on FHIR
- Smart Backend Services App (Multi-patient authorization and API) : These are server-to-server backend applications e.g. Export clinical data for multiple patients in a group. This app is system level app without any UI.

Types of Authentications supported

- Symmetric ("client secret") authentication
- Asymmetric ("private key JWT") authentication
- Public Clients

Steps for Smart App Launch

- Client App Registration
- Retrieve .well-known/smart-configuration
- Obtain authorization code
- Obtain access token
- Access FHIR API
- Refresh Token to obtain new access token

Client APP Registration

Before a SMART app can run against an EHR, the app must be registered with that EHR's authorization service. We are using OAuth 2.0 Client Registration.

- Client App Request process
- Third-party application must sign the API Subscription Agreement with following details
 - The app name
 - The necessary APIs / Scopes
 - Any redirect URIs
 - Launch_URL (optional)
 - URL to JWK Set (Only for apps supporting asymmetric client authentication).
- The application is created by the EHR admin after confirming the app's registration parameters and communicates a client_id to the app.
 - FHIR APIs (of any supported version) listed within the USCDI v3 core data set will be supported
 - Only reads data from FHIR server
- Retrieve .well-known/smart-configuration

In order to obtain launch context and request authorization to access FHIR resources, the app discovers the EHR FHIR server's SMART configuration metadata, including OAuth authorization_endpoint and token_endpoint URLs.

FHIR server makes SMART configuration available from well-known endpoint. You can get Authorization end point and token endpoint.

Request

GET 200 https://fhir.iPatientCare.net:9443/fhirserver/fhir/.well-known/smart-configuration

Response

```
{  
  "capabilities": [  

```

```

"launch-ehr",
"launch-standalone",
"client-public",
"client-confidential-symmetric",
"sso-openid-connect",
"context-banner",
"context-style",
"context-ehr-patient",
"context-ehr-encounter",
"context-standalone-patient",
"context-standalone-encounter",
"permission-offline",
"permission-patient",
"permission-user"
],
"revocation_endpoint": "https://oauth.iPatientCare.net/outh/token/revoke-token",
"authorization_endpoint": "https://oauth.iPatientCare.net/outh/authorization",
"token_endpoint": "https://oauth.iPatientCare.net/outh/token"
}

```

Error codes

JSON response is returned from the FHIR Server along with the code

200: Ok

401: Unauthorized

400: Bad Request

500: Internal Server Error

Obtain authorization code

The app supplies the following parameters to the EHR's "authorize" endpoint.

Parameters

`response_ty` `require` Fixed value: `code`.

Parameters

<code>client_id</code>	require	The client's identifier.
<code>redirect_uri</code>	require	Must match one of the client's pre-registered redirect URIs.
<code>launch</code>	condition	When using the EHR Launch flow, this must match the launch value received from the EHR. Omitted when using the Standalone Launch .
<code>scope</code>	require	Must describe the access that the app needs, including scopes like <code>patient/*.rs</code>, <code>openid</code> and <code>fhirUser</code> (if app needs authenticated patient identity) and either: - a <code>launch</code> value indicating that the app wants to receive already-established launch context details from the EHR - a set of launch context requirements in the form <code>launch/patient</code>, which asks the EHR to establish context on your behalf. Please refer to the scopes supported in the table below
<code>state</code>	require	An opaque value used by the client to maintain state between the request and callback. The authorization server includes this value when redirecting the user-agent back to the client. The parameter SHALL be used for preventing cross-site request forgery or session fixation attacks. The app SHALL use an unpredictable value for the state parameter with at least 122 bits of entropy (e.g., a properly configured random uuid is suitable).
<code>aud</code>	require	URL of the EHR resource server from which the app wishes to retrieve FHIR data. This parameter prevents leaking a genuine bearer token to a counterfeit resource server. (Note that in the case of an EHR launch flow, this <code>aud</code> value is the same as the launch's <code>iss</code> value.) Note that the <code>aud</code> parameter is semantically equivalent to the <code>resource</code> parameter defined in RFC8707 . SMART's <code>aud</code> parameter predates RFC8707 and we have decided not to rename it for reasons of backwards compatibility. We may consider renaming SMART's <code>aud</code> parameter in the future if implementer feedback indicates that alignment would be valuable. For the current release, servers SHALL support the <code>aud</code> parameter and MAY support a <code>resource</code> parameter as a synonym for

Scopes Supported

Scope	Grants
patient/*.i	Permission to read and search any resource for the current patient (see notes on wildcard scopes below).
user/*.rea	Permission to read and write all resources that the current user can access (see notes on wildcard scopes below).
openid fhirUser	Permission to retrieve information about the current logged-in user.
launch	Permission to obtain launch context when app is launched from an EHR.
launch/pa	When launching outside the EHR, ask for a patient to be selected at launch time.
offline_ac	Request a refresh_token that can be used to obtain a new access token to replace an expired one, even after the end-user no longer is online after the access token expires.

Here's an example of an authorization request using HTTP GET. You will replace the [redirect_uri], [client_id], [launch_token], [state], [code_challenge], and [audience] placeholders with your own values.

Request

https://ehr/authorize?

response_type=code&
 client_id=app-client-id&
 redirect_uri=https%3A%2F%2Fapp%2Fafter-auth&
 launch=xyz123&
 scope=launch+patient%2FObservation.rs+patient%2FPatient.rs+openid+fhirUser&
 state=98wrghuwuogerg97&
 aud=https://ehr/fhir

e.g.

https://oauth.iPatientCare.net/outh/authorization?response_type=code&client_id=SAMPLE_CONFIDENTIAL_CLIENT_ID&redirect_uri=https%3A%2F%2Finferno.healthit.gov%2Finferno%2Foauth2%2Fstatic%2Fredirect&scope=launch%2Fpatient+openid+fhirUser+offline_access+patient%2FMedication.read+patient%2FAllergyIntolerance.read+patient%2FCarePlan.read+patient%2FCareTeam.read+patient%2FCondition.read+patient%2FDevice.read+patient%2FDiagnosticReport.read+patient%2FDocumentReference.read+patient%2FEncounter.read+patient%2FGoal.read+patient%2FImmunization.read+patient%2FLocation.read+patient%2FMedicationRequest.read+patient%2FObservation.read+patient%2FOrganization.read+patient%2FPatient.r

ead+patient%2FPractitioner.read+patient%2FProcedure.read+patient%2FProvenance.read+patient%2FPractionerRole.read&state=26a45e37-7445-4e4a-b8fb-24144eccbdc4&aud=https%3A%2F%2Finferno.healthit.gov%2Freference-server%2Fr4

Response

The EHR's authorization server reviews the request from your application. If approved, the authorization server redirects the browser to the redirect URL supplied in the initial request and appends the following querystring parameter.

- **code:** This parameter contains the authorization code generated by EHR, which will be exchanged for the access token in the next step.
- **state:** The exact value received from the client.

Here is how redirect URL will look like

Location: <https://app/after-auth?code=123abc&state=98wrghuwoogerg97>

https://inferno.healthit.gov/inferno/oauth2/static/redirect?code=SAMPLE_CODE.bGF1bmNoL3BhdGllbnQgb3BlbmlkIGZoaXJvc2VyIG9mZmxpbmVfYWNjZXNzIHdhbGllbnQvTWVkaWNhdGlvbi5yZWZkIHdhbGllbnQvQWxsZXJneUludG9sZXJhbmNlLnJlYWQgcGF0aWVudC9DYXJlUGxhbi5yZWZkIHdhbGllbnQvQ2FyZVRIYW0ucmVhZCBwYXRpZW50L0NvbRpdGlvbi5yZWZkIHdhbGllbnQvRGV2aWNlLnJlYWQgcGF0aWVudC9EaWFnbm9zdGljUmVwb3J0LnJlYWQgcGF0aWVudC9Eb2N1bWVudFJlZmVyZW5jZS5yZWZkIHdhbGllbnQvRW5jb3Vu dGVyLnJlYWQgcGF0aWVudC9Hb2FsLnJlYWQgcGF0aWVudC9JbW11bml6YXRpb24ucmVhZCBwYXRpZW50 L0xvY2F0aW9uLnJlYWQgcGF0aWVudC9NZWRpY2F0aW9uUmVxdWVzdC5yZWZkIHdhbGllbnQvT2JzZXJ2YX Rpb24ucmVhZCBwYXRpZW50L09yZ2FuaXphdGlvbi5yZWZkIHdhbGllbnQvUGF0aWVudC5yZWZkIHdhbGllbn QvUHJhY3RpdGlvbmVyLnJlYWQgcGF0aWVudC9Qcm9jZWR1cmUucmVhZCBwYXRpZW50L1Byb3ZlbnFuY2 UucmVhZCBwYXRpZW50L1ByYWN0aXRpb25lclJvbGUucmVhZCA=.ODU=&state=26a45e37-7445-4e4a-b8fb-24144eccbdc4

After receiving the authorization code, your application trades the code for a JSON object containing an access token and contextual information by sending an HTTP POST to the token endpoint using a Content-Type header with value of "application/x-www-form-urlencoded".

Error codes

JSON response is returned from the FHIR Server along with the code

200: Ok

401: Unauthorized

400: Bad Request

500: Internal Server Error

Obtain access token

For public apps, authentication not required because a client with no secret cannot prove its identity when it issues a call. (The end-to-end system can still be secure because the client comes from a known, https protected endpoint specified and enforced by the redirect uri.) For confidential apps, authentication is required. Confidential clients SHOULD use Asymmetric Authentication if available, and MAY use Symmetric Authentication.

There are two different requests to get access token based on authentication type. In case of Symmetric ("client secret") authentication app issues an HTTP POST to the EHR authorization server's token endpoint URL using content-type application/x-www-form-urlencoded.

In case of Asymmetric authentication app, generate a client authentication assertion and prepare arguments for POST to token API:

Request

Parameters		
grant_type	required	Fixed value: authorization_code for symmetric and client_credentials for assymetric
Code	required	Code that the app received from the authorization server
redirect_uri	required	The same redirect_uri used in the initial authorization request
client_id	condition	Required for public apps . Omit for confidential apps .
client_assertion_t	condition	Required for asymmetric authentication. set to urn:ietf:params:oauth:client-assertion-type:jwt-bearer
client_assertion	condition	Required for asymmetric authentication. set to a JWT signed with dynamic client's private key

Parameters		
Scope	condition	system/*.read for backend services type of application (bulk operations)

Symmetric Authentication Request

<https://oauth.iPatientCare.net/oauth/token>

Payload

grant_type=authorization_code&code=SAMPLE_CODE.bGF1bmNoL3BhdGllbnQgb3BlbmklGZoaXJvc2VylG9mZmxpbmVfYWNjZXNzIHBhdGllbnQvTWVkaWNhdGlvbi5yZWFKIHBhdGllbnQvQWxsZXJneUludG9sZXJhbmlnLnJlYWQgcGF0aWVudC9DYXJlUGxhbi5yZWFKIHBhdGllbnQvQ2FyZVRIYW0ucmVhZCBwYXRpZW50L0NvbWdRpdGlvi5yZWFKIHBhdGllbnQvRGV2aWNLnJlYWQgcGF0aWVudC9EaWFnbm9zdGljUmVwb3J0LnJlYWQgcGF0aWVudC9Eb2N1bWVudFJlZmVyZW5jZS5yZWFKIHBhdGllbnQvRW5jb3VudGVyLnJlYWQgcGF0aWVudC9Hb2FsLnJlYWQgcGF0aWVudC9JbW11bml6YXRpb24ucmVhZCBwYXRpZW50L0xvY2F0aW9uLnJlYWQgcGF0aWVudC9NZWRpY2F0aW9uUmVxdWVzdC5yZWFKIHBhdGllbnQvT2JzZXJ2YXRpb24ucmVhZCBwYXRpZW50L09yZ2FuaXphdGlvi5yZWFKIHBhdGllbnQvUGF0aWVudC5yZWFKIHBhdGllbnQvUHJhY3RpdGlvmVylLnJlYWQgcGF0aWVudC9Qcm9jZWR1cmUucmVhZCBwYXRpZW50L1Byb3ZlbnFuY2UucmVhZCBwYXRpZW50L1ByYWN0aXRpb25lclJvbGUucmVhZCA%3D.ODU%3D&redirect_uri=https%3A%2F%2Finferno.healthit.gov%2Finferno%2Foauth%2Fstatic%2Fredirect

In case of Asymmetric authentication following is a request

client_assertion=eyJ0eXAiOiJKV1QiLCJhbGciOiJIUzI1NiIsImtpZCI6IjRiNDIhNzU5ZDFiYjExNWlzMjI1ZjRjZjliZWl2ZDFlbn0.eyJpc3MiOiJleUpoYkdjaU9pSkIvekkxTmIjc0luUjVjQ0k2SWtwWFZDSXNjbXRwWkNjNkluSmxammx6ZUhKaGRHbHZiaTEwYjJ0bGpSjkuZXIKcWQydHpYM1Z5YkNjNkltadBkSEE2THk4eE1DNHhOUzR5TIRjdU56TXZhVzVtWlHkdWJ5OHVkbMlZzYkMxcmlJOTNiaTlxZDJ0ekxtcHpiMjRpbTENKaFkyTmxjM05VYjJ0bGJuTkZISEJwY21WSmJpSTZNVFVzSW1saGRDSTZNVFU1TnpReE16RTVOWDAucTR2NE1zYzc0a041MDZLVFowcV9taW55YXBkdzBnd2xUNk1fdWIMNzNTNClInN1Yil6ImV5SmhiR2NpT2IKSVV6STFOaUlzSW5SNWNDSTZJa3BYVkJ0c0ltdHBaQ0k2SW5KbFoybHpkSEpoZEdsdmJpMTBiMnRsYmlKOS5leUpxZDJ0elgzVnliQ0k2SW1oMGRIQTZMeTh4TUM0eE5TNHIOVEI1TnpNdmFXNW1aWEp1Ynk4dWQyVnNiQzFyYm05M2JpOXFkMnR6TG1wemlyNGIMQ0poWTJObGMzTlViMnRsYm5ORmVIQnBjbVZKYmlJNk1UVXNjbWxoZENjNk1UVTV0eIF4TXpFNU5YMC5xNHY0TXNjNzRrTjUwNktUWjBxX21pbnlhcEp3MGd3bFQ2TV91aUw3M1M0liwiYXVkljoiaHR0cHM6Ly9pbmZlcm5vLmhlYWx0aGl0Lmdvdi9yZWZlcmVuY2Utc2VydmVyl29hdXR0L2J1bGstdG9rZW4iLCJleHAiOiJlE2NTMxMDY4ODgsImp0aSI6ImZmZjA5YjlxMjhiOWRmNmU3MGU5MWQ2ZDUwMGQxMmFmYmU2YWE4N2FhMjJlZjVhMjc5YjZhNTI0YzA0NmM1ZWEifQ.tlx7jlu9wn7QdHHdWyUMf9otY0r7eX62LuVF3kZOsjiw-CP-glqdbDtOKiu-BQ9doduW4zyk3fs6k48XpAEQLsSXajfEJeDy9uj-WdSmpAHEC0xYHQsZajakdWDSBRA&client_assertion_type=urn%3Aietf%3Aparams%3Aoauth%3Aclient-assertion-type%3Ajwt-bearer&grant_type=client_credentials&scope=system%2F%2A.read

Response

The EHR authorization server SHALL return a JSON object that includes an access token or a message indicating that the authorization request has been denied. The JSON structure includes the following parameters:

Parameters

access_token	required	The access token issued by the authorization server
token_type	required	Fixed value: Bearer
expires_in	recommended	Lifetime in seconds of the access token, after which the token SHALL NOT be accepted by the resource server
scope	required	Scope of access authorized. Note that this can be different from the scopes requested by the app.
id_token	optional	Authenticated user identity and user details, if requested
refresh_token	optional	Token that can be used to obtain a new access token, using the same or a subset of the original authorization grants

Sample response

[illegible]

ounE6oLyYZJp_OGvVhkREu8jO77m3Tgsji_jbX7g4-
deuwB4F9EHUfpjfhM3TB1GloQ7cgJFhHOs9mTCcpT0aYYQkmCJGRx1R5jl-I56p4_63IkU2BRau-
INII3Zvcnsz5ajpvU46elwg",

"smart_style_url": "https://oauth.iPatientCare.net/app/smart-style-url",

"token_type": "bearer",

"expires_in": 3600

}

Error codes

JSON response is returned from the FHIR Server along with the code

200: Ok

401: Unauthorized

400: Bad Request

500: Internal Server Error

At this point, the authorization flow is complete.

Access FHIR API

With a valid access token, the app can access protected EHR data by issuing a FHIR API call to the FHIR endpoint on the EHR's resource server. Only HTTP GET request types are supported while accessing FHIR Resource API. Following resources are supported

FHIR Resource

Allergies and Intolerances

Assessment and Plan of Treatment

Care Team Members

Clinical Notes -US Core DocumentReference Profile&US Core DiagnosticReport Profile for Report and I Exchange

Goals - US Core Goal Profile

Health Concerns - US Core Condition Profile

Immunizations - US Core Immunization Profile

Laboratory - US Core Laboratory Result Observation Profile & US Core DiagnosticReport Profile for Laboratory Results Reporting

Medications - US Core Medication Profile, US Core MedicationRequest Profile, US Core AllergyIntolerance Profile

Patient Demographics - US Core Patient Profile

Problems - US Core Condition Profile

Procedures - US Core Procedure Profile

Provenance - US Core Provenance Profile

Smoking Status - US Core Smoking Status Observation Profile

Unique Device Identifier(s) for Patient's Implantable Device(s) - US Core ImplantableDevice Profile

Vital Signs—FHIR Core VitalSigns Profile

Organization – US Core Organization Profile

DocumentReference – US Core DocumentReference Profile

ServiceRequest – US Core ServiceRequest Profile

Coverage – US Core Coverage Profile

MedicationDispense – US Core MedicationDispense Profile

Specimen – US Core Specimen Profile

RelatedPerson – US Core RelatedPerson Profile

To make the FHIR API call to the FHIR endpoint, you need to pass the Tenant-id of the practice in the request header "X-FHIR-TENANT-ID". Tenant-id is a unique code assigned by AssureCare to each practice on the FHIR Server. Please contact the practice to obtain their Tenant-id.

Request : Patient

Search Parameters

_id: GET [Base_URL]/Patient/_search?Patient=[id]

Identifier : GET [Base_URL]/Patient?identifier={system}}[code]

Name: GET [Base_URL]/Patient?name=[string]

birthdate + name :GET [Base_URL]/Patient?birthdate=[date]&name=[string]

gender + name : GET [Base_URL]/Patient?gender={system}}[code]&name=[string]

Response

US Core Patient Profile as specified in the [US Core v3.1.1 Implementation Guide](#)

Code: 200 (Success)

```
{
  "resourceType": "Bundle",
  "id": "6f1d6f2d-c42b-453b-8f5a-8d74c24f6fb6",
  "type": "searchset",
  "total": 1,
  "link": [
    {
      "relation": "self",
      "url": "https://FHIR-Base-URL/Patient?_count=10&_id=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&_page=1"
    }
  ],
  "entry": [
    {
      "id": "17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b",
```

```

"fullUrl": "https://FHIR-Base-URL/Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b",
"resource": {
  "resourceType": "Patient",
  "id": "17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b",
  "meta": {
    "versionId": "3",
    "lastUpdated": "2022-03-24T10:50:49.231Z",
    "profile": [
      "http://hl7.org/fhir/us/core/StructureDefinition/us-core-patient"
    ]
  },
  "text": {
    "status": "generated",
    "div": "<div xmlns=\"http://www.w3.org/1999/xhtml\">Generated by <a href=\"https://github.com/synthetichealth/synthea\">Synthea</a>.Version identifier: v2.6.1-3-g50f4f58f\n . Person seed: 4286870567281389426 Population seed: 3</div>"
  },
  "extension": [
    {
      "extension": [
        {
          "url": "ombCategory",
          "valueCoding": {
            "system": "urn:oid:2.16.840.1.113883.6.238",
            "code": "2106-3",
            "display": "White"
          }
        },
        {
          "url": "text",
          "valueString": "White"
        }
      ]
    },
    {
      "url": "http://hl7.org/fhir/us/core/StructureDefinition/us-core-race"
    }
  ],

```

```

{
  "extension": [
    {
      "url": "ombCategory",
      "valueCoding": {
        "system": "urn:oid:2.16.840.1.113883.6.238",
        "code": "2186-5",
        "display": "Non Hispanic or Latino"
      }
    },
    {
      "url": "text",
      "valueString": "Non Hispanic or Latino"
    }
  ],
  "url": "http://hl7.org/fhir/us/core/StructureDefinition/us-core-ethnicity",
},
{
  "url": "http://hl7.org/fhir/StructureDefinition/patient-mothersMaidenName",
  "valueString": "Christen366 Murray856"
},
{
  "url": "http://hl7.org/fhir/us/core/StructureDefinition/us-core-birthsex",
  "valueCode": "M"
},
{
  "url": "http://hl7.org/fhir/StructureDefinition/patient-birthPlace",
  "valueAddress": {
    "city": "Springfield",
    "state": "Massachusetts",
    "country": "US"
  }
}
],

```

```

"identifier": [
  {
    "system": "https://github.com/synthetichealth/synthea",
    "value": "e91975f5-9445-c11f-cabf-c3c6dae161f2"
  },
  {
    "type": {
      "coding": [
        {
          "system": "http://terminology.hl7.org/CodeSystem/v2-0203",
          "code": "MR",
          "display": "Medical Record Number"
        }
      ],
      "text": "Medical Record Number"
    },
    "system": "http://hospital.smarthealthit.org",
    "value": "e91975f5-9445-c11f-cabf-c3c6dae161f2"
  },
  {
    "type": {
      "coding": [
        {
          "system": "http://terminology.hl7.org/CodeSystem/v2-0203",
          "code": "SS",
          "display": "Social Security Number"
        }
      ],
      "text": "Social Security Number"
    },
    "system": "http://hl7.org/fhir/sid/us-ssn",
    "value": "999-61-9798"
  },
  {

```

```

"type": {
  "coding": [
    {
      "system": "http://terminology.hl7.org/CodeSystem/v2-0203",
      "code": "DL",
      "display": "Driver's License"
    }
  ],
  "text": "Driver's License"
},
"system": "urn:oid:2.16.840.1.113883.4.3.25",
"value": "S99940282"
}
],
"name": [
  {
    "use": "official",
    "family": "Harry",
    "given": [
      "Tom"
    ],
    "prefix": [
      "Mr."
    ],
    "suffix": [
      "Jr."
    ]
  },
  {
    "use": "maiden",
    "given": [
      "Alex"
    ]
  }
]
}

```

```

],
"telecom": [
  {
    "system": "phone",
    "value": "555-770-2787",
    "use": "home"
  }
],
"gender": "male",
"birthDate": "1941-09-05",
"address": [
  {
    "extension": [
      {
        "extension": [
          {
            "url": "latitude",
            "valueDecimal": 42.1343534042923
          },
          {
            "url": "longitude",
            "valueDecimal": -72.67217549422628
          }
        ]
      },
      {
        "url": "http://hl7.org/fhir/StructureDefinition/geolocation"
      }
    ]
  },
  {
    "line": [
      "599 Schowalter Promenade"
    ],
    "city": "West Springfield",
    "state": "MA",
    "postalCode": "01089",
    "country": "US",

```

```

    "period": {
      "start": "1940-09-05"
    }
  },
  "maritalStatus": {
    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/v3-MaritalStatus",
        "code": "M",
        "display": "M"
      }
    ],
    "text": "M"
  },
  "multipleBirthBoolean": false,
  "communication": [
    {
      "language": {
        "coding": [
          {
            "system": "urn:ietf:bcp:47",
            "code": "en-US",
            "display": "English"
          }
        ],
        "text": "English"
      }
    }
  ],
  "search": {
    "mode": "match",
    "score": 1
  }
}

```

```
}  
}  
]  
}
```

Request : AllergyIntolerance

Search Parameters

GET GET [Base_URL]/AllergyIntolerance?Patient=[id]

Response

Each resource returned from the first search is expected to conform to the [US Core AllergyIntolerance Profile](#).

Code: 200 (Success)

```
{
  "resourceType": "AllergyIntolerance",
  "id": "17e7758c319-7955ecfe-94da-4071-8488-4ef21425ed85",
  "meta": {
    "versionId": "1",
    "lastUpdated": "2022-01-20T11:56:52.889Z",
    "profile": [
      "http://hl7.org/fhir/us/core/StructureDefinition/us-core-allergyintolerance"
    ]
  },
  "clinicalStatus": {
    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/allergyintolerance-clinical",
        "code": "inactive"
      }
    ]
  },
  "verificationStatus": {
    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/allergyintolerance-verification",
        "code": "confirmed"
      }
    ]
  }
}
```

```

},
"type": "allergy",
"category": [
  "food"
],
"criticality": "low",
"code": {
  "coding": [
    {
      "system": "http://snomed.info/sct",
      "code": "300916003",
      "display": "Latex allergy"
    }
  ],
  "text": "Latex allergy"
},
"patient": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdfb2be4b"
},
"recordedDate": "1941-10-14T18:33:18-05:00",
"reaction": [
  {
    "manifestation": [
      {
        "coding": [
          {
            "system": "http://snomed.info/sct",
            "code": "271807003",
            "display": "skin rash"
          }
        ],
        "text": "skin rash"
      }
    ]
  }
]

```

```
}  
]  
}
```

Request : CarePlan

GET [Base_URL]/AllergyIntolerance?CarePlan=assess-plan&patient=[id]

e.g.

[https://\[Base_URL\]/CarePlan?category=assess-plan&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b](https://[Base_URL]/CarePlan?category=assess-plan&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b)

Response

Each resource returned from the first search is expected to conform to the [US Core CarePlan Profile](#).

Code: 200 (Success)

```
{
  "resourceType": "Bundle",
  "id": "79d7a78e-de9b-4c5c-b51d-45b72963ae1f",
  "type": "searchset",
  "total": 1,
  "link": [
    {
      "relation": "self",
      "url": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/CarePlan?_count=10&category=assess-plan&patient=Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&_page=1"
    }
  ],
  "entry": [
    {
      "id": "17e776e13ff-1ea7505d-f475-4f49-adde-30e69df86734",
      "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/CarePlan/17e776e13ff-1ea7505d-f475-4f49-adde-30e69df86734",
      "resource": {
        "resourceType": "CarePlan",
        "id": "17e776e13ff-1ea7505d-f475-4f49-adde-30e69df86734",
        "meta": {
          "versionId": "2",
          "lastUpdated": "2022-01-21T10:56:06.304Z",
```

```

    "profile": [
      "http://hl7.org/fhir/us/core/StructureDefinition/us-core-careplan"
    ]
  },
  "text": {
    "status": "generated",
    "div": "<div xmlns=\"http://www.w3.org/1999/xhtml\">Care Plan for Fracture care.<br/>Activities:
<ul><li>Fracture care</li><li>Fracture care</li></ul><br/>Care plan is meant to treat Fracture of
clavicle.</div>"
  },
  "status": "active",
  "intent": "order",
  "category": [
    {
      "coding": [
        {
          "system": "http://hl7.org/fhir/us/core/CodeSystem/careplan-category",
          "code": "assess-plan"
        }
      ]
    },
    {
      "coding": [
        {
          "system": "http://snomed.info/sct",
          "code": "384758001",
          "display": "Self-care interventions (procedure)"
        }
      ],
      "text": "Self-care interventions (procedure)"
    }
  ],
  "subject": {
    "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
  },

```

```

"encounter": {
  "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
},
"period": {
  "start": "1941-10-02T18:33:18-05:00"
},
"careTeam": [
  {
    "reference": "CareTeam/17e775e7546-36e6a4cd-264e-4a94-92fa-cdf1b2fa30c1"
  }
],
"activity": [
  {
    "detail": {
      "code": {
        "coding": [
          {
            "system": "http://snomed.info/sct",
            "code": "409002",
            "display": "Food allergy diet"
          }
        ],
        "text": "Food allergy diet"
      },
      "status": "in-progress",
      "location": {
        "reference": "Location/17e7750e68f-6266aec8-0cd3-443e-adb9-be37f3bf25f2",
        "display": "HOLYOKE MEDICAL CENTER"
      }
    }
  },
  {
    "detail": {
      "code": {

```

```

    "coding": [
      {
        "system": "http://snomed.info/sct",
        "code": "58332002",
        "display": "Allergy education"
      }
    ],
    "text": "Allergy education"
  },
  "status": "in-progress",
  "location": {
    "reference": "Location/17e7750e68f-6266aec8-0cd3-443e-adb9-be37f3bf25f2",
    "display": "HOLYOKE MEDICAL CENTER"
  }
},
{
  "detail": {
    "code": {
      "coding": [
        {
          "system": "http://snomed.info/sct",
          "code": "58332002",
          "display": "Allergy education"
        }
      ],
      "text": "Allergy education"
    },
    "status": "in-progress",
    "location": {
      "reference": "Location/17e7750e68f-6266aec8-0cd3-443e-adb9-be37f3bf25f2",
      "display": "HOLYOKE MEDICAL CENTER"
    }
  }
}

```

```
    }  
  ]  
},  
"search": {  
  "mode": "match",  
  "score": 1  
}  
}  
]  
}
```

Request : CareTeam

GET [Base_URL]/CareTeam?patient=[id]&status=active

e.g.

https://[Base_URL]/CareTeam?patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&status=active

Response

Each resource returned from the first search is expected to conform to the [US Core CareTeam Profile](#)

Code: 200 (Success)

```
{
  "resourceType": "Bundle",
  "id": "492509fd-0650-4e2e-834f-41fc7b132a5e",
  "type": "searchset",
  "total": 1,
  "link": [
    {
      "relation": "self",
      "url": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/CareTeam?_count=10&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&status=http://hl7.org/fhir/care-team-status%7Cactive&_page=1"
    }
  ],
  "entry": [
    {
      "id": "17e775e7546-36e6a4cd-264e-4a94-92fa-cdf1b2fa30c1",
      "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/CareTeam/17e775e7546-36e6a4cd-264e-4a94-92fa-cdf1b2fa30c1",
      "resource": {
        "resourceType": "CareTeam",
        "id": "17e775e7546-36e6a4cd-264e-4a94-92fa-cdf1b2fa30c1",
        "meta": {
          "versionId": "1",
```

```

    "lastUpdated": "2022-01-20T12:03:06.182Z",
    "profile": [
      "http://hl7.org/fhir/us/core/StructureDefinition/us-core-careteam"
    ]
  },
  "status": "active",
  "subject": {
    "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
  },
  "encounter": {
    "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
  },
  "period": {
    "start": "1941-10-02T18:33:18-05:00"
  },
  "participant": [
    {
      "role": [
        {
          "coding": [
            {
              "system": "http://snomed.info/sct",
              "code": "116154003",
              "display": "Patient"
            }
          ],
          "text": "Patient"
        }
      ],
      "member": {
        "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b",
        "display": "Mr. Dustin31 Ritchie586"
      }
    }
  ],

```

```

{
  "role": [
    {
      "coding": [
        {
          "system": "http://snomed.info/sct",
          "code": "223366009",
          "display": "Healthcare professional (occupation)"
        }
      ],
      "text": "Healthcare professional (occupation)"
    }
  ],
  "member": {
    "reference": "Practitioner/17e774bb850-6b7b7f50-ed52-4523-a3a8-60fb0ef71c5b",
    "display": "Dr. Patricia625 Kilback373"
  }
},
{
  "role": [
    {
      "coding": [
        {
          "system": "http://snomed.info/sct",
          "code": "224891009",
          "display": "Healthcare services (qualifier value)"
        }
      ],
      "text": "Healthcare services (qualifier value)"
    }
  ],
  "member": {
    "reference": "Organization/17e774dca91-e7e4e9e3-390e-4c10-83ae-88beabc176e2",
    "display": "HOLYOKE MEDICAL CENTER"
  }
}

```

```
    }
  }
],
"managingOrganization": [
  {
    "reference": "Organization/17e774dca91-e7e4e9e3-390e-4c10-83ae-88beabc176e2",
    "display": "HOLYOKE MEDICAL CENTER"
  }
]
},
"search": {
  "mode": "match",
  "score": 1
}
}
]
```

Request: Problems / Health Concern

Search Parameters

GET [base]/Condition?patient=[id]&category={system|}[code]

GET [base]/Condition?patient=[id]&category=http://terminology.hl7.org/CodeSystem/condition-category|problem-list-item

GET [base]/Condition?patient=[id]&category=http://hl7.org/fhir/us/core/CodeSystem/condition-category|health-concern

GET [base]/Condition?patient=[id]&category=http://terminology.hl7.org/CodeSystem/condition-category|encounter-diagnosis

e.g

[Base_URL]/Condition?patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&category=http://hl7.org/fhir/us/core/CodeSystem/condition-category|health-concern

Response

Each resource returned from the first search is expected to conform to the US Core [Condition](#) Profile

Code: 200 (Success)

```
{
  "resourceType": "Bundle",
  "id": "f537c30e-fab7-45f0-8167-4ab7a48386ec",
  "type": "searchset",
  "total": 2,
  "link": [
    {
      "relation": "self",
      "url": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Condition?_count=10&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&_page=1"
    }
  ],
  "entry": [
    {
      "id": "17e7777ce8c-c66752fb-ee1b-40c2-ba5a-2e6ee89b3e1a",
      "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Condition/17e7777ce8c-c66752fb-ee1b-40c2-ba5a-2e6ee89b3e1a",

```

```

"resource": {
  "resourceType": "Condition",
  "id": "17e7777ce8c-c66752fb-ee1b-40c2-ba5a-2e6ee89b3e1a",
  "meta": {
    "versionId": "1",
    "lastUpdated": "2022-01-20T12:30:47.436Z",
    "profile": [
      "http://hl7.org/fhir/us/core/StructureDefinition/us-core-condition"
    ]
  },
  "clinicalStatus": {
    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/condition-clinical",
        "code": "resolved"
      }
    ]
  },
  "verificationStatus": {
    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/condition-ver-status",
        "code": "confirmed"
      }
    ]
  },
  "category": [
    {
      "coding": [
        {
          "system": "http://terminology.hl7.org/CodeSystem/condition-category",
          "code": "encounter-diagnosis",
          "display": "Encounter Diagnosis"
        }
      ]
    }
  ]
}

```

```

    ]
  }
],
"code": {
  "coding": [
    {
      "system": "http://snomed.info/sct",
      "code": "444814009",
      "display": "Viral sinusitis (disorder)"
    }
  ],
  "text": "Viral sinusitis (disorder)"
},
"subject": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
},
"encounter": {
  "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
},
"onsetDateTime": "1940-11-09T18:33:18-05:00",
"abatementDateTime": "1940-11-30T18:33:18-05:00",
"recordedDate": "1940-11-09T18:33:18-05:00"
},
"search": {
  "mode": "match",
  "score": 1
}
},
{
  "id": "17eb029f925-40603327-0ac3-48dd-972b-6c882b8751f5",
  "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Condition/17eb029f925-40603327-0ac3-48dd-972b-6c882b8751f5",
  "resource": {
    "resourceType": "Condition",
    "id": "17eb029f925-40603327-0ac3-48dd-972b-6c882b8751f5",

```

```

"meta": {
  "versionId": "1",
  "lastUpdated": "2022-01-31T12:44:07.847Z",
  "profile": [
    "http://hl7.org/fhir/us/core/StructureDefinition/us-core-condition"
  ]
},
"clinicalStatus": {
  "coding": [
    {
      "system": "http://terminology.hl7.org/CodeSystem/condition-clinical",
      "code": "resolved"
    }
  ]
},
"verificationStatus": {
  "coding": [
    {
      "system": "http://terminology.hl7.org/CodeSystem/condition-ver-status",
      "code": "confirmed"
    }
  ]
},
"category": [
  {
    "extension": [
      {
        "url": "http://hl7.org/fhir/StructureDefinition/data-absent-reason",
        "valueCode": "unknown"
      }
    ]
  }
],
"code": {

```

```

    "coding": [
      {
        "system": "http://snomed.info/sct",
        "code": "195662009",
        "display": "Acute viral pharyngitis (disorder)"
      }
    ],
    "text": "Acute viral pharyngitis (disorder)"
  },
  "subject": {
    "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
  },
  "encounter": {
    "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
  },
  "onsetDateTime": "1961-06-25T01:11:45-04:00",
  "abatementDateTime": "1961-07-07T01:11:45-04:00",
  "recordedDate": "1961-06-25T01:11:45-04:00"
},
"search": {
  "mode": "match",
  "score": 1
}
]
}

```

Request: Implantable Device

Search Parameters

GET [Base_URL]/Device?patient=[id]

e.g. [Base_URL]/Device?patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b

Response

Each resource returned from the first search is expected to conform to the US Core [ImplantableDevice](#) Profile

Code: 200 (Success)

```
{
  "resourceType": "Bundle",
  "id": "5804551e-120d-4164-b05b-9705fad44fa5",
  "type": "searchset",
  "total": 1,
  "link": [
    {
      "relation": "self",
      "url": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Device?_count=10&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&_page=1"
    }
  ],
  "entry": [
    {
      "id": "17e777b3f34-04bc3627-92d3-43bb-92e7-35533110f750",
      "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Device/17e777b3f34-04bc3627-92d3-43bb-92e7-35533110f750",
      "resource": {
        "resourceType": "Device",
        "id": "17e777b3f34-04bc3627-92d3-43bb-92e7-35533110f750",
        "meta": {
          "versionId": "1",
          "lastUpdated": "2022-01-20T12:34:32.884Z",
```

```

"profile": [
  "http://hl7.org/fhir/us/core/StructureDefinition/us-core-implantable-device"
],
"udiCarrier": [
  {
    "deviceIdentifier": "43069338026389",
    "carrierHRF": "(01)43069338026389(11)000302(17)250317(10)1134(21)842026117977"
  }
],
"status": "active",
"distinctIdentifier": "43069338026389",
"manufactureDate": "2000-03-02T18:33:18-05:00",
"expirationDate": "2025-03-17T19:33:18-04:00",
"lotNumber": "1134",
"serialNumber": "842026117977",
"deviceName": [
  {
    "name": "Implantable defibrillator, device (physical object)",
    "type": "user-friendly-name"
  }
],
"type": {
  "coding": [
    {
      "system": "http://snomed.info/sct",
      "code": "72506001",
      "display": "Implantable defibrillator, device (physical object)"
    }
  ],
  "text": "Implantable defibrillator, device (physical object)"
},
"patient": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
}

```

```
    }  
  },  
  "search": {  
    "mode": "match",  
    "score": 1  
  }  
}  
]  
}
```

Request: Diagnostic Report/ Clinical Notes/ DocumentReference

Search Parameters

patient : GET [Base_URL]/DiagnosticReport?patient=[id]

patient + code : GET [Base_URL]/DiagnosticReport?patient=[id]&code={system}}[code]

patient + category: GET [Base_URL]/DiagnosticReport?patient=[id]&category=LAB

patient : GET [Base_URL]/ DocumentReference?patient=[id]

e.g.

[https://\[Base_URL\]/DiagnosticReport?category=LP29684-5&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b](https://[Base_URL]/DiagnosticReport?category=LP29684-5&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b)

[https://\[Base_URL\]/DiagnosticReport?category=LP7839-6&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b](https://[Base_URL]/DiagnosticReport?category=LP7839-6&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b)

[https://\[Base_URL\]/DiagnosticReport?category=LP29708-2&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b](https://[Base_URL]/DiagnosticReport?category=LP29708-2&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b)

[https://\[Base_URL\]/DocumentReference?patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b](https://[Base_URL]/DocumentReference?patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b)

Category and codes details can be obtained from [US Core DiagnosticReport Profile for Report and Note exchange](#)

Response

Each resource returned from the first search is expected to conform to the [US Core DiagnosticReport Profile for Report and Note exchange](#).

Code: 200 (Success)

```
{
  "resourceType": "Bundle",
  "id": "a87a8318-2dfc-488d-9f0a-6b07d7c29271",
  "type": "searchset",
  "total": 1,
  "link": [
    {
      "relation": "self",
```

```

    "url":
    "https://fhir.iPatientCare.net:9443/fhirserver/fhir/DiagnosticReport?_count=10&category=LP29684-
    5&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&_page=1"
  }
],
"entry": [
  {
    "id": "17e7b6e6f1c-119e9b8b-11b9-448b-b760-dc116013b654",
    "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/DiagnosticReport/17e7b6e6f1c-119e9b8b-
    11b9-448b-b760-dc116013b654",
    "resource": {
      "resourceType": "DiagnosticReport",
      "id": "17e7b6e6f1c-119e9b8b-11b9-448b-b760-dc116013b654",
      "meta": {
        "versionId": "1",
        "lastUpdated": "2022-01-21T06:59:02.053Z",
        "profile": [
          "http://hl7.org/fhir/us/core/StructureDefinition/us-core-diagnosticreport-note"
        ]
      },
      "status": "final",
      "category": [
        {
          "coding": [
            {
              "system": "http://loinc.org",
              "code": "LP29684-5",
              "display": "Radiology"
            }
          ],
          "text": "Radiology"
        }
      ],
      "code": {
        "coding": [

```

```

{
  "system": "http://loinc.org",
  "code": "34117-2",
  "display": "History and physical note"
}
],
"text": "History and physical note"
},
"subject": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
},
"encounter": {
  "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
},
"effectiveDateTime": "1974-01-24T19:33:18-04:00",
"issued": "1974-01-24T19:33:18.715-04:00",
"performer": [
  {
    "reference": "Practitioner/17e774bb850-6b7b7f50-ed52-4523-a3a8-60fb0ef71c5b",
    "display": "Dr. Blossom971 Christiansen251"
  },
  {
    "reference": "Practitioner/17e774bb850-6b7b7f50-ed52-4523-a3a8-60fb0ef71c5b"
  },
  {
    "reference": "Organization/17e774dca91-e7e4e9e3-390e-4c10-83ae-88beabc176e2"
  }
],
"presentedForm": [
  {
    "contentType": "text/plain",
    "data":
"CJE5NzQtMDEtMjQKCiMgQ2hpZWYgQ29tcGxhaW50Ci0gRGVjcmVhc2VklGlulEp1ZGdlbWVudAotlENvbm
Z1c2lvbG9tIE5hc2FsIEVnbmdlc3Rpb24KLSBTbmVlemluZyBGaXRzCi0gTmFzYWwgRGZlY2hhcmdlCi0gQ291Z
2gKCgojIEhpc3Rvcnkgb2YgUHJlc2VudCBJbGxuZXNzCkR1c3RpbjMxIGZlIGegMzMgeWVhcn1vbGQgbm9uL
Whpc3BhbmljIHdoaXRlIG1hbGUuIFBhdGllbnQgaGFzIGEgaGlzdG9yeSBvZiBhY3V0ZSBhbGxlcmdpYyByZWJj

```

dGlvbiwgYWN1dGUgYnJvbmNoaXRpcyAoZGlzb3JkZXlpLCBsYWNlcmF0aW9uIG9mIHROaWdoLCBvdGl0aX
MgbWVkaWEsIHZpcmFslHNpbnVzaXRpcyAoZGlzb3JkZXlpLCBzcHJhaW4gb2YgYW5rbGUslGZpcnN0IGRIZ3
JlZSBidXJuLCBmcmFjdHVyZSBvZiBhbmtsZSwgY2hpbGRob29klGFzdGhtYS4KCIMgU29jaWFsIEhpc3RvcnkKU
GF0aWVudCBpcyBtYXJyaWVvLiBQYXRpZW50IGlzlHRlbXBvcnFyaWx5IGhvbWVsZXNzLiBQYXRpZW50IGlzlG
FulGFjdGl2ZSBzbW9rZXIgaW5kIGlzlGFuIGFsY29ob2xpYy4gUGF0aWVudCBpZGVudGlmaWVzIGFzIGhldGVy
b3NleHVhbC4KCIBhdGllbnQgY29tZXMGZnJvbSBhIG1pZGRsZSBzb2Npb2Vjb25vbWljIGJhY2tncm91bmQuIF
BhdGllbnQgZGlkIG5vdCBmaW5pc2ggaGlnaCBzY2hvb2wuIFBhdGllbnQgY3VycmVudGx5IGhhcyBBZXRUYS4
KCIMgQWxsZXJnaWVzCmxhdGV4IGFsbGVyZ3ksIGFsbGVyZ3kgdG8gc295YQoKlyBNZWRpY2F0aW9ucwphY
2V0YW1pbm9waGVuIDMyNSBtZyBvcnFslHRhYmxldDsgMTlwIGFjdHVhdCBmbHV0aWNhc29uZSBwcm9wa
W9uYXRlIDAuMDQ0IG1nL2FjdHVhdCBtZXRIcmVklGRvc2UgaW5oYWxlcjsgbmRhMDlwNTAzIDlwMCBhY3R1
YXQgYWxidXRlcm9sIDAuMDkgaWVvYWN0dWF0IG1ldGVyZWQgZG9zZSBpbmhhbGVyOyBwcmVkbmlzb25l
IDUgbWcgY3JhbCB0YWJsZXQ7IGFzcGlyaW4gODEgbWcgY3JhbCB0YWJsZXQKCIMgQXNzZXNzbWVudCBh
bmQgUGxhbgoKCIMjIFBsYW4KUGF0aWVudCB3YXMgZ2l2ZW4gdGhlIGZvbGxvd2luZyBpbW11bml6YXRpb2
5zOiBpbmZsdWVuemEslHNIYXNvbmFsLCBpbmplY3RhYmxlLCBwcmVzZXJ2YXRpdmlUgZnJlZS4gCIRoZSBm
b2xs3dpbmcmcgHJvY2VkdXJlcyB3ZXJlIGNvbmlR1Y3RIZDoKLSBtZWRpY2F0aW9uIHJlY29uY2lsaWF0aW9uCh
wcm9jZWR1cmUpCIRoZSBwYXRpZW50IHdhcyBwcmVzY3JpYmVklHRoZSBmb2xs3dpbmcmgbWVkaWNhdG
lvbnM6Ci0gMTlwIGFjdHVhdCBmbHV0aWNhc29uZSBwcm9waW9uYXRlIDAuMDQ0IG1nL2FjdHVhdCBtZXRI
cmVklGRvc2UgaW5oYWxlcgotIG5kYTAyMDUwMyAyMDAgYWN0dWF0IGFsYnV0ZXJvbCAwLjA5IG1nL2FjdH
VhdCBtZXRIcmVklGRvc2UgaW5oYWxlcgo="

```
    }  
  ]  
},  
"search": {  
  "mode": "match",  
  "score": 1  
}  
}  
]  
}
```

Request: Laboratory Result Observation

Search Parameters

patient : GET [Base_URL]/DiagnosticReport?patient=[id]

patient + code : GET [Base_URL]/DiagnosticReport?patient=[id]&code={system}[code]

patient + category: GET [Base_URL]/DiagnosticReport?patient=[id]&category=LAB

e.g.

`https://[Base_URL]/DiagnosticReport?category=http://terminology.hl7.org/CodeSystem/v2-0074%7CLAB&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b`

`https://[Base_URL]/DiagnosticReport?category=LAB&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b`

`https://[Base_URL]/DiagnosticReport?category=LAB&patient=Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b`

Category and codes details can be obtained from [US Core DiagnosticReport Profile for Laboratory Results Reporting](#).

Response

Each resource returned from the first search is expected to conform to the [US Core DiagnosticReport Profile for Laboratory Results Reporting](#)

Code: 200 (Success)

```
{
  "resourceType": "Bundle",
  "id": "019290b1-c4e1-4025-862c-60ffc3a6160b",
  "type": "searchset",
  "total": 1,
  "link": [
    {
      "relation": "self",
      "url":
        "https://fhir.iPatientCare.net:9443/fhirserver/fhir/DiagnosticReport?_count=10&category=http://terminology.hl7.org/CodeSystem/v2-0074%7CLAB&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&_page=1"
    }
  ],
}
```

```

"entry": [
  {
    "id": "17e7b7b7e66-36f4cf84-3482-4e36-8cca-13cc832b3bed",
    "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/DiagnosticReport/17e7b7b7e66-36f4cf84-3482-4e36-8cca-13cc832b3bed",
    "resource": {
      "resourceType": "DiagnosticReport",
      "id": "17e7b7b7e66-36f4cf84-3482-4e36-8cca-13cc832b3bed",
      "meta": {
        "versionId": "1",
        "lastUpdated": "2022-01-21T07:13:17.926Z",
        "profile": [
          "http://hl7.org/fhir/us/core/StructureDefinition/us-core-diagnosticreport-lab"
        ]
      },
      "status": "final",
      "category": [
        {
          "coding": [
            {
              "system": "http://terminology.hl7.org/CodeSystem/v2-0074",
              "code": "LAB",
              "display": "Laboratory"
            }
          ]
        }
      ],
      "code": {
        "coding": [
          {
            "system": "http://loinc.org",
            "code": "51990-0",
            "display": "Basic Metabolic Panel"
          }
        ]
      }
    }
  ],

```

```

    "text": "Basic Metabolic Panel"
  },
  "subject": {
    "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
  },
  "encounter": {
    "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
  },
  "effectiveDateTime": "2017-09-28T19:33:18-04:00",
  "issued": "2017-09-28T19:33:18.715-04:00",
  "performer": [
    {
      "reference": "Organization/17e774dca91-e7e4e9e3-390e-4c10-83ae-88beabc176e2",
      "display": "PCP170967"
    },
    {
      "reference": "Practitioner/17e774bb850-6b7b7f50-ed52-4523-a3a8-60fb0ef71c5b"
    },
    {
      "reference": "Organization/17e774dca91-e7e4e9e3-390e-4c10-83ae-88beabc176e2"
    }
  ],
  "result": [
    {
      "reference": "Observation/17e77809e90-d5cfd26d-ab5d-4f36-bc32-1ee799b2c278",
      "display": "Glucose"
    },
    {
      "reference": "Observation/17e7785ea69-6c9f88d9-0553-4216-be17-f908ebe0b371",
      "display": "Urea Nitrogen"
    }
  ]
},
"search": {

```

```
    "mode": "match",  
    "score": 1  
  }  
}  
]  
}
```

Request: Goal

GET https://[Base_URL]/Goal?patient=[id]

e.g. https://[Base_URL]/Goal?patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b

Response

Each resource returned from the first search is expected to conform to the [US Core Goal Profile](#)

Code: 200 (Success)

```
{
  "resourceType": "Bundle",
  "id": "6f2db230-b4f4-4447-856b-3f398377a95f",
  "type": "searchset",
  "total": 1,
  "link": [
    {
      "relation": "self",
      "url": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Goal?_count=10&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&_page=1"
    }
  ],
  "entry": [
    {
      "id": "17e7b8b62f8-a4f77df6-b97e-4ece-82ff-36480c9cc38c",
      "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Goal/17e7b8b62f8-a4f77df6-b97e-4ece-82ff-36480c9cc38c",
      "resource": {
        "resourceType": "Goal",
        "id": "17e7b8b62f8-a4f77df6-b97e-4ece-82ff-36480c9cc38c",
        "meta": {
          "versionId": "1",
          "lastUpdated": "2022-01-21T07:30:39.48Z",
          "profile": [
            "http://hl7.org/fhir/us/core/StructureDefinition/us-core-goal"
          ]
        }
      }
    }
  ]
}
```

```

},
"lifecycleStatus": "active",
"description": {
  "coding": [
    {
      "system": "http://snomed.info/sct",
      "code": "281004",
      "display": "Alcoholic dementia"
    }
  ],
  "text": "Alcoholic dementia"
},
"subject": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdfb2be4b"
},
"target": [
  {
    "dueDate": "2020-12-07"
  }
]
},
"search": {
  "mode": "match",
  "score": 1
}
}
]
}

```

Request: Immunization

GET https://[Base_URL]/Immunization?patient=[id]

e.g. https://[Base_URL]/Immunization?patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b

Response

Each resource returned from the first search is expected to conform to the US Core [Immunization](#) Profile
Code: 200 (Success)

```
{
  "resourceType": "Bundle",
  "id": "206f31dc-9636-4e6f-a652-40c307eebe40",
  "type": "searchset",
  "total": 2,
  "link": [
    {
      "relation": "self",
      "url":
"https://fhir.iPatientCare.net:9443/fhirserver/fhir/Immunization?_count=10&patient=Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&_page=1"
    }
  ],
  "entry": [
    {
      "id": "17e7b8dc9ab-8abc85e3-0ab1-470f-ae8-1f4e9a3597bb",
      "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Immunization/17e7b8dc9ab-8abc85e3-0ab1-470f-ae8-1f4e9a3597bb",
      "resource": {
        "resourceType": "Immunization",
        "id": "17e7b8dc9ab-8abc85e3-0ab1-470f-ae8-1f4e9a3597bb",
        "meta": {
          "versionId": "1",
          "lastUpdated": "2022-01-21T07:33:16.843Z",
          "profile": [
            "http://hl7.org/fhir/us/core/StructureDefinition/us-core-immunization"
          ]
        },

```

```

"status": "not-done",
"statusReason": {
  "coding": [
    {
      "system": "http://terminology.hl7.org/CodeSystem/v3-ActReason",
      "code": "OSTOCK",
      "display": "product out of stock"
    }
  ],
  "text": "product out of stock"
},
"vaccineCode": {
  "coding": [
    {
      "system": "http://terminology.hl7.org/CodeSystem/data-absent-reason",
      "code": "unknown",
      "display": "Unknown"
    }
  ],
  "text": "Unknown"
},
"patient": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
},
"encounter": {
  "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
},
"occurrenceString": "1946-08-29T19:33:18-04:00",
"primarySource": true,
"location": {
  "reference": "Location/17e7750e68f-6266aec8-0cd3-443e-adb9-be37f3bf25f2",
  "display": "RIVER BEND MEDICAL GROUP - CHICOPEE OFFICE - URGENT CARE"
},
},

```

```

"search": {
  "mode": "match",
  "score": 1
},
{
  "id": "17eb02c0f47-3b8a6c05-ae47-4002-9001-fd61626bfca5",
  "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Immunization/17eb02c0f47-3b8a6c05-ae47-4002-9001-fd61626bfca5",
  "resource": {
    "resourceType": "Immunization",
    "id": "17eb02c0f47-3b8a6c05-ae47-4002-9001-fd61626bfca5",
    "meta": {
      "versionId": "1",
      "lastUpdated": "2022-01-31T12:46:24.583Z",
      "profile": [
        "http://hl7.org/fhir/us/core/StructureDefinition/us-core-immunization"
      ]
    },
    "status": "not-done",
    "statusReason": {
      "coding": [
        {
          "system": "http://terminology.hl7.org/CodeSystem/v3-ActReason",
          "code": "OSTOCK",
          "display": "product out of stock"
        }
      ],
      "text": "product out of stock"
    },
    "vaccineCode": {
      "coding": [
        {
          "system": "http://terminology.hl7.org/CodeSystem/data-absent-reason",
          "code": "unknown",

```

```

      "display": "Unknown"
    }
  ],
  "text": "Unknown"
},
"patient": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
},
"encounter": {
  "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
},
"occurrenceString": "1946-08-29T19:33:18-04:00",
"primarySource": true,
"location": {
  "reference": "Location/17e7750e68f-6266aec8-0cd3-443e-adb9-be37f3bf25f2",
  "display": "RIVER BEND MEDICAL GROUP - CHICOPEE OFFICE - URGENT CARE"
}
},
"search": {
  "mode": "match",
  "score": 1
}
}
]
}

```

Request: Medication

GET [https://\[Base_URL\]/MedicationRequest?patient=\[id\]](https://[Base_URL]/MedicationRequest?patient=[id])

e.g. [https://\[Base_URL\]/MedicationRequest?patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b &intent=order&status=active](https://[Base_URL]/MedicationRequest?patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&intent=order&status=active)

Response

Each resource returned from the first search is expected to conform to the US Core [Medication](#) Profile, US Core [MedicationRequest](#) Profile

Code: 200 (Success)

```
{
  "resourceType": "MedicationRequest",
  "id": "17e7b9b263e-9cd8f64f-f7e1-447e-8583-a753342d13c6",
  "meta": {
    "versionId": "1",
    "lastUpdated": "2022-01-21T07:47:52.51Z",
    "profile": [
      "http://hl7.org/fhir/us/core/StructureDefinition/us-core-medicationrequest"
    ]
  },
  "status": "stopped",
  "intent": "plan",
  "reportedBoolean": true,
  "medicationCodeableConcept": {
    "coding": [
      {
        "system": "http://www.nlm.nih.gov/research/umls/rxnorm",
        "code": "2123111",
        "display": "NDA020503 200 ACTUAT Albuterol 0.09 MG/ACTUAT Metered Dose Inhaler"
      }
    ],
    "text": "NDA020503 200 ACTUAT Albuterol 0.09 MG/ACTUAT Metered Dose Inhaler"
  },
  "subject": {
    "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
  },
}
```

```

"authoredOn": "2020-08-13T19:33:18-04:00",
"requester": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
},
"reasonReference": [
  {
    "reference": "Condition/17e7777ce8c-c66752fb-ee1b-40c2-ba5a-2e6ee89b3e1a"
  }
],
"dosageInstruction": [
  {
    "sequence": 1,
    "text": "Take as needed.",
    "asNeededBoolean": true
  }
]
}

```

Request: Smoking Status

GET [https://\[Base_URL\]/Observation?code=72166-2&patient=\[id\]](https://[Base_URL]/Observation?code=72166-2&patient=[id])

e.g. [https://\[Base_URL\]/Observation?code=72166-2&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b](https://[Base_URL]/Observation?code=72166-2&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b)

Response

Each resource returned from the first search is expected to conform to the US Core [Smoking Status Observation](#) Profile

Code: 200 (Success)

```

{
  "resourceType": "Bundle",
  "id": "2ea85043-b450-40f1-8d2d-751af1c10abb",
  "type": "searchset",
  "total": 2,
  "link": [
    {

```

```

    "relation": "self",
    "url": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Observation?_count=10&code=72166-2&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&_page=1"
  }
],
"entry": [
  {
    "id": "17e77809e90-d5cfd26d-ab5d-4f36-bc32-1ee799b2c278",
    "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Observation/17e77809e90-d5cfd26d-ab5d-4f36-bc32-1ee799b2c278",
    "resource": {
      "resourceType": "Observation",
      "id": "17e77809e90-d5cfd26d-ab5d-4f36-bc32-1ee799b2c278",
      "meta": {
        "versionId": "1",
        "lastUpdated": "2022-01-20T12:40:24.976Z",
        "profile": [
          "http://hl7.org/fhir/us/core/StructureDefinition/us-core-smokingstatus"
        ]
      },
      "status": "final",
      "category": [
        {
          "coding": [
            {
              "system": "http://terminology.hl7.org/CodeSystem/observation-category",
              "code": "survey",
              "display": "survey"
            }
          ]
        }
      ],
      "code": {
        "coding": [
          {

```

```

    "system": "http://loinc.org",
    "code": "72166-2",
    "display": "Tobacco smoking status NHIS"
  }
],
"text": "Tobacco smoking status NHIS"
},
"subject": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
},
"encounter": {
  "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
},
"effectivePeriod": {
  "start": "1940-09-05T19:33:18-04:00",
  "end": "1940-09-05T20:33:18-04:00"
},
"issued": "1940-09-05T19:33:18.715-04:00",
"valueCodeableConcept": {
  "coding": [
    {
      "system": "http://snomed.info/sct",
      "code": "266919005",
      "display": "Never smoker"
    }
  ],
  "text": "Never smoker"
}
},
"search": {
  "mode": "match",
  "score": 1
}
},

```

```

{
  "id": "17e7785ea69-6c9f88d9-0553-4216-be17-f908ebe0b371",
  "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Observation/17e7785ea69-6c9f88d9-0553-4216-be17-f908ebe0b371",
  "resource": {
    "resourceType": "Observation",
    "id": "17e7785ea69-6c9f88d9-0553-4216-be17-f908ebe0b371",
    "meta": {
      "versionId": "1",
      "lastUpdated": "2022-01-20T12:46:12.073Z",
      "profile": [
        "http://hl7.org/fhir/us/core/StructureDefinition/us-core-smokingstatus"
      ]
    },
    "status": "final",
    "category": [
      {
        "coding": [
          {
            "system": "http://terminology.hl7.org/CodeSystem/observation-category",
            "code": "survey",
            "display": "survey"
          }
        ]
      }
    ],
    "code": {
      "coding": [
        {
          "system": "http://loinc.org",
          "code": "72166-2",
          "display": "Tobacco smoking status NHIS"
        }
      ],
      "text": "Tobacco smoking status NHIS"
    }
  }
}

```

```

},
"subject": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882dfdb2be4b"
},
"encounter": {
  "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
},
"effectivePeriod": {
  "start": "1940-09-05T19:33:18-04:00",
  "end": "1940-09-05T20:33:18-04:00"
},
"issued": "1940-09-05T19:33:18.715-04:00",
"valueCodeableConcept": {
  "coding": [
    {
      "system": "http://snomed.info/sct",
      "code": "266919005",
      "display": "Never smoker"
    }
  ],
  "text": "Never smoker"
}
},
"search": {
  "mode": "match",
  "score": 1
}
}
]
}

```

Request: Procedure

GET https://[Base_URL]/Procedure?patient=[id]

e.g. https://[Base_URL]/Procedure?patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b

Response

Each resource returned from the first search is expected to conform to the US Core [Procedure](#) Profile

Code: 200 (Success)

```
{
  "resourceType": "Bundle",
  "id": "ed4425dd-7479-4b52-b404-502697281ee0",
  "type": "searchset",
  "total": 1,
  "link": [
    {
      "relation": "self",
      "url": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Procedure?_count=10&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&_page=1"
    }
  ],
  "entry": [
    {
      "id": "17e7ba24dee-bf67ef6b-6419-4c7c-b905-a107b26bf8fb",
      "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Procedure/17e7ba24dee-bf67ef6b-6419-4c7c-b905-a107b26bf8fb",
      "resource": {
        "resourceType": "Procedure",
        "id": "17e7ba24dee-bf67ef6b-6419-4c7c-b905-a107b26bf8fb",
        "meta": {
          "versionId": "1",
          "lastUpdated": "2022-01-21T07:55:41.422Z",
          "profile": [
            "http://hl7.org/fhir/us/core/StructureDefinition/us-core-procedure"
          ]
        },
        "status": "completed",
```

```

"code": {
  "coding": [
    {
      "system": "http://snomed.info/sct",
      "code": "430193006",
      "display": "Medication Reconciliation (procedure)"
    }
  ],
  "text": "Medication Reconciliation (procedure)"
},
"subject": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
},
"encounter": {
  "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
},
"performedDateTime": "1941-02-13T18:33:18-05:00",
"location": {
  "reference": "Location/17e7750e68f-6266aec8-0cd3-443e-adb9-be37f3bf25f2",
  "display": "PCP170967"
}
},
"search": {
  "mode": "match",
  "score": 1
}
]
}

```

Request: Provenance

GET https://[Base_URL]/Provenance/[id]

e.g.

https://[Base_URL]/Provenance/17e8adbdd12-21675728-3f52-4320-8b3b-d5977915336c

Response

Each resource returned from the first search is expected to conform to the US Core [Provenance](#) Profile
Code: 200 (Success)

```
{
  "resourceType": "Provenance",
  "id": "17e8adbdd12-21675728-3f52-4320-8b3b-d5977915336c",
  "meta": {
    "versionId": "5",
    "lastUpdated": "2022-01-31T12:51:12.934Z",
    "profile": [
      "http://hl7.org/fhir/us/core/StructureDefinition/us-core-provenance"
    ]
  },
  "target": [
    {
      "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdfb2be4b"
    },
    {
      "reference": "Location/17e7750e68f-6266aec8-0cd3-443e-adb9-be37f3bf25f2"
    },
    {
      "reference": "Organization/17e774dca91-e7e4e9e3-390e-4c10-83ae-88beabc176e2"
    },
    {
      "reference": "Practitioner/17e774bb850-6b7b7f50-ed52-4523-a3a8-60fb0ef71c5b"
    },
    {
      "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
    }
  ]
}
```

```

{
  "reference": "AllergyIntolerance/17e7758c319-7955ecfe-94da-4071-8488-4ef21425ed85"
},
{
  "reference": "CareTeam/17e775e7546-36e6a4cd-264e-4a94-92fa-cdf1b2fa30c1"
},
{
  "reference": "CarePlan/17e776e13ff-1ea7505d-f475-4f49-adde-30e69df86734"
},
{
  "reference": "Condition/17e7777ce8c-c66752fb-ee1b-40c2-ba5a-2e6ee89b3e1a"
},
{
  "reference": "Condition/17eb029f925-40603327-0ac3-48dd-972b-6c882b8751f5"
},
{
  "reference": "Device/17e777b3f34-04bc3627-92d3-43bb-92e7-35533110f750"
},
{
  "reference": "Observation/17e77809e90-d5cfd26d-ab5d-4f36-bc32-1ee799b2c278"
},
{
  "reference": "Observation/17e7785ea69-6c9f88d9-0553-4216-be17-f908ebe0b371"
},
{
  "reference": "Observation/17e778a9113-66c19e50-1b3f-41fd-a4da-ef86393fb0b2"
},
{
  "reference": "Observation/17e778d42db-92528d6f-7103-40fb-99ca-b2ac8837210d"
},
{
  "reference": "Observation/17e7791dd19-11eaf5c1-a5c0-429c-9867-f3110c743572"
},
{

```

```

    "reference": "Observation/17e77947034-7677a2e6-33af-4e06-83ae-a6727ac79885"
  },
  {
    "reference": "Observation/17e7797685b-e707227d-0cb5-4b28-bab2-ca1f4a1ba111"
  },
  {
    "reference": "Observation/17e7799a564-8f9ed437-d396-46bf-8f68-bf728955654c"
  },
  {
    "reference": "Observation/17e779bc672-37c76705-1f73-4531-bfcb-ee3cd3e25bfd"
  },
  {
    "reference": "Observation/17e77a065f3-77947ea0-1b53-47a5-9a9e-10f6e169f292"
  },
  {
    "reference": "Observation/17e77a29d41-85ec6814-a99a-4a73-ae02-f6a494a11468"
  },
  {
    "reference": "Observation/17e77a53c30-b356fe0d-d90b-42ed-91f9-5ed08c8047ca"
  },
  {
    "reference": "Observation/17e7c606d8a-9aeab44c-0261-465d-8930-4a218b201bf5"
  },
  {
    "reference": "DiagnosticReport/17e7b6e6f1c-119e9b8b-11b9-448b-b760-dc116013b654"
  },
  {
    "reference": "DiagnosticReport/17e7b741f29-4b6b813d-d11e-4597-8178-5bc5996ac4e4"
  },
  {
    "reference": "DiagnosticReport/17e7b7b7e66-36f4cf84-3482-4e36-8cca-13cc832b3bed"
  },
  {
    "reference": "DiagnosticReport/17e7d0e6d1d-c0e0136b-81d4-4164-a8b8-863694276a4f"
  }

```

```

},
{
  "reference": "DocumentReference/17e7b838c61-4be99b80-0f94-4b26-a086-a30ed78c91ce"
},
{
  "reference": "DocumentReference/17e7d144d28-7d79ed3e-ca9e-400d-a9ed-add23f5cef07"
},
{
  "reference": "DocumentReference/17e7d19866c-dbbaa2e9-4d81-44bb-af5d-ec2e829e9c72"
},
{
  "reference": "DocumentReference/17e7d1c49cd-bfb30b25-9f98-40c7-9d94-977b7b803f21"
},
{
  "reference": "DocumentReference/17e7d24c947-205cd8e5-3a93-4ab9-8662-5e92b0080af7"
},
{
  "reference": "Goal/17e7b8b62f8-a4f77df6-b97e-4ece-82ff-36480c9cc38c"
},
{
  "reference": "Immunization/17e7b8dc9ab-8abc85e3-0ab1-470f-ae8-1f4e9a3597bb"
},
{
  "reference": "Immunization/17eb02c0f47-3b8a6c05-ae47-4002-9001-fd61626bfca5"
},
{
  "reference": "Medication/17e7b8f47d0-2bf2943a-123e-4725-b5dd-f5200df6c3a2"
},
{
  "reference": "MedicationRequest/17e7b94dd8d-d3b0e43a-de55-4822-8bb5-07e5fd6177d9"
},
{
  "reference": "MedicationRequest/17e7b990c43-4525f6a4-81a0-40f7-97c7-7dbebff6e34b"
},

```

```

{
  "reference": "MedicationRequest/17e7b9b263e-9cd8f64f-f7e1-447e-8583-a753342d13c6"
},
{
  "reference": "MedicationRequest/17e7ba05cb1-0a5defbd-deb9-4792-9952-8f7f0e5bc74b"
},
{
  "reference": "Procedure/17e7ba24dee-bf67ef6b-6419-4c7c-b905-a107b26bf8fb"
},
{
  "reference": "Observation/17e81e97c4b-ea761488-0b67-43d0-a0d7-346dea549bee"
},
{
  "reference": "Observation/17e81ebac04-75d9f0f2-fcb2-4e2f-9c22-87fb7ddf5f5e"
},
{
  "reference": "Observation/17e81ed9c56-bcda0427-df92-4b58-b3ae-ac11c22f4d7d"
},
{
  "reference": "Observation/17e81ef5b5e-40df3e7d-19bc-4c2a-94bc-9a04ff93537c"
},
{
  "reference": "Observation/17e81f0fd34-6fc95b32-8346-4aec-9de3-072d12c992a5"
},
{
  "reference": "Observation/17e81f291e7-c99c91ac-57dd-47f6-b7f7-a365d2c087d7"
},
{
  "reference": "Observation/17e81f4a61d-ff3926f2-e35b-466a-bb6c-62a061425258"
},
{
  "reference": "Observation/17e8501692e-8c2d145f-fb69-4ae3-85d1-ef31ebf4a6a2"
},
{

```

```

    "reference": "Observation/17e8502f180-3c7e4770-3630-4654-922a-3f802c2939b1"
  }
],
"recorded": "2020-12-10T18:33:18.715-05:00",
"agent": [
  {
    "type": {
      "coding": [
        {
          "system": "http://terminology.hl7.org/CodeSystem/provenance-participant-type",
          "code": "author",
          "display": "Author"
        }
      ],
      "text": "Author"
    },
    "who": {
      "reference": "Practitioner/17e774bb850-6b7b7f50-ed52-4523-a3a8-60fb0ef71c5b",
      "display": "Dr. Blossom971 Christiansen251"
    },
    "onBehalfOf": {
      "reference": "Organization/17e774dca91-e7e4e9e3-390e-4c10-83ae-88beabc176e2",
      "display": "PCP170967"
    }
  },
  {
    "type": {
      "coding": [
        {
          "system": "http://hl7.org/fhir/us/core/CodeSystem/us-core-provenance-participant-type",
          "code": "transmitter",
          "display": "Transmitter"
        }
      ],

```

```
    "text": "Transmitter"
  },
  "who": {
    "reference": "Practitioner/17e774bb850-6b7b7f50-ed52-4523-a3a8-60fb0ef71c5b",
    "display": "Dr. Blossom971 Christiansen251"
  },
  "onBehalfOf": {
    "reference": "Organization/17e774dca91-e7e4e9e3-390e-4c10-83ae-88beabc176e2",
    "display": "PCP170967"
  }
}
]
```

Request: Vital Signs

patient + code: GET [https://\[Base_URL\]/Observation?code=85354-9&patient=\[id\]](https://[Base_URL]/Observation?code=85354-9&patient=[id])

patient + category + date : GET [https://\[Base_URL\]/Observation?category=vital-signs&date=\[leYYYY-MM-DD\]&patient=\[id\]](https://[Base_URL]/Observation?category=vital-signs&date=[leYYYY-MM-DD]&patient=[id])

e.g.

[https://\[Base_URL\]/Observation?code=85354-9&patient=Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b](https://[Base_URL]/Observation?code=85354-9&patient=Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b)

[https://\[Base_URL\]/Observation?code=85354-9&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b](https://[Base_URL]/Observation?code=85354-9&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b)

[https://\[Base_URL\]/Observation?category=vital-signs&date=gt1994-05-18T19:33:18-04:00&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b](https://[Base_URL]/Observation?category=vital-signs&date=gt1994-05-18T19:33:18-04:00&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b)

Response

Each resource returned from the first search is expected to conform to the FHIR Core [VitalSigns](#) Profile

Code: 200 (Success)

```
{
  "resourceType": "Bundle",
  "id": "231df7ad-8e46-4371-99cc-4ce0dc9c5cb0",
  "type": "searchset",
  "total": 6,
  "link": [
    {
      "relation": "self",
      "url": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Observation?_count=10&category=vital-signs&date=gt1994-05-18T19:33:18-04:00&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&_page=1"
    }
  ],
  "entry": [
    {
      "id": "17e7797685b-e707227d-0cb5-4b28-bab2-ca1f4a1ba111",
      "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Observation/17e7797685b-e707227d-0cb5-4b28-bab2-ca1f4a1ba111",
      "resource": {
        "resourceType": "Observation",
```

```

"id": "17e7797685b-e707227d-0cb5-4b28-bab2-ca1f4a1ba111",
"meta": {
  "versionId": "1",
  "lastUpdated": "2022-01-20T13:05:18.427Z",
  "profile": [
    "http://hl7.org/fhir/StructureDefinition/bodytemp",
    "http://hl7.org/fhir/StructureDefinition/vitalsigns"
  ]
},
"status": "final",
"category": [
  {
    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/observation-category",
        "code": "vital-signs",
        "display": "vital-signs"
      }
    ]
  }
],
"code": {
  "coding": [
    {
      "system": "http://loinc.org",
      "code": "8310-5",
      "display": "Body temperature"
    },
    {
      "system": "http://loinc.org",
      "code": "8331-1",
      "display": "Oral temperature"
    }
  ],

```

```

    "text": "Body temperature"
  },
  "subject": {
    "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
  },
  "encounter": {
    "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
  },
  "effectiveDateTime": "2009-08-19T19:33:18-04:00",
  "issued": "2009-08-19T19:33:18.715-04:00",
  "valueQuantity": {
    "value": 37.799,
    "unit": "Cel",
    "system": "http://unitsofmeasure.org",
    "code": "Cel"
  }
},
"search": {
  "mode": "match",
  "score": 1
}
},
{
  "id": "17e7799a564-8f9ed437-d396-46bf-8f68-bf728955654c",
  "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Observation/17e7799a564-8f9ed437-d396-46bf-8f68-bf728955654c",
  "resource": {
    "resourceType": "Observation",
    "id": "17e7799a564-8f9ed437-d396-46bf-8f68-bf728955654c",
    "meta": {
      "versionId": "1",
      "lastUpdated": "2022-01-20T13:07:45.124Z",
      "profile": [
        "http://hl7.org/fhir/StructureDefinition/bp",
        "http://hl7.org/fhir/StructureDefinition/vitalsigns"
      ]
    }
  }
}

```

```

    ]
  },
  "status": "final",
  "category": [
    {
      "coding": [
        {
          "system": "http://terminology.hl7.org/CodeSystem/observation-category",
          "code": "vital-signs",
          "display": "vital-signs"
        }
      ]
    }
  ],
  "code": {
    "coding": [
      {
        "system": "http://loinc.org",
        "code": "85354-9",
        "display": "Blood Pressure"
      }
    ],
    "text": "Blood Pressure"
  },
  "subject": {
    "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
  },
  "encounter": {
    "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
  },
  "effectiveDateTime": "1994-05-19T19:33:18-04:00",
  "issued": "1994-05-19T19:33:18.715-04:00",
  "component": [
    {

```

```

"code": {
  "coding": [
    {
      "system": "http://loinc.org",
      "code": "8462-4",
      "display": "Diastolic Blood Pressure"
    }
  ],
  "text": "Diastolic Blood Pressure"
},
"valueQuantity": {
  "value": 85,
  "unit": "mm[Hg]",
  "system": "http://unitsofmeasure.org",
  "code": "mm[Hg]"
}
},
{
  "code": {
    "coding": [
      {
        "system": "http://loinc.org",
        "code": "8480-6",
        "display": "Systolic Blood Pressure"
      }
    ],
    "text": "Systolic Blood Pressure"
  },
  "valueQuantity": {
    "value": 129,
    "unit": "mm[Hg]",
    "system": "http://unitsofmeasure.org",
    "code": "mm[Hg]"
  }
}

```

```

    }
  ]
},
"search": {
  "mode": "match",
  "score": 1
}
},
{
  "id": "17e77a29d41-85ec6814-a99a-4a73-ae02-f6a494a11468",
  "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Observation/17e77a29d41-85ec6814-a99a-4a73-ae02-f6a494a11468",
  "resource": {
    "resourceType": "Observation",
    "id": "17e77a29d41-85ec6814-a99a-4a73-ae02-f6a494a11468",
    "meta": {
      "versionId": "1",
      "lastUpdated": "2022-01-20T13:17:32.865Z",
      "profile": [
        "http://hl7.org/fhir/StructureDefinition/hearttrate",
        "http://hl7.org/fhir/StructureDefinition/vitalsigns"
      ]
    },
    "status": "final",
    "category": [
      {
        "coding": [
          {
            "system": "http://terminology.hl7.org/CodeSystem/observation-category",
            "code": "vital-signs",
            "display": "vital-signs"
          }
        ]
      }
    ]
  },

```

```

"code": {
  "coding": [
    {
      "system": "http://loinc.org",
      "code": "8867-4",
      "display": "Heart rate"
    }
  ],
  "text": "Heart rate"
},
"subject": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882dfdb2be4b"
},
"encounter": {
  "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
},
"effectiveDateTime": "2003-07-10T19:33:18-04:00",
"issued": "2003-07-10T19:33:18.715-04:00",
"valueQuantity": {
  "value": 65,
  "unit": "/min",
  "system": "http://unitsofmeasure.org",
  "code": "/min"
},
},
"search": {
  "mode": "match",
  "score": 1
}
},
{
  "id": "17e81ef5b5e-40df3e7d-19bc-4c2a-94bc-9a04ff93537c",
  "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Observation/17e81ef5b5e-40df3e7d-19bc-4c2a-94bc-9a04ff93537c",
  "resource": {

```

```

"resourceType": "Observation",
"id": "17e81ef5b5e-40df3e7d-19bc-4c2a-94bc-9a04ff93537c",
"meta": {
  "versionId": "1",
  "lastUpdated": "2022-01-22T13:17:34.43Z",
  "source": "#iLFSV7OLv0KF8dmQ",
  "profile": [
    "http://hl7.org/fhir/StructureDefinition/vitalsigns",
    "http://hl7.org/fhir/StructureDefinition/bodytemp"
  ]
},
"status": "final",
"category": [
  {
    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/observation-category",
        "code": "vital-signs",
        "display": "vital-signs"
      }
    ]
  }
],
"code": {
  "coding": [
    {
      "system": "http://loinc.org",
      "code": "8310-5",
      "display": "Body temperature"
    },
    {
      "system": "http://loinc.org",
      "code": "8331-1",
      "display": "Oral temperature"
    }
  ]
}

```

```

    }
  ],
  "text": "Body temperature"
},
"subject": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
},
"encounter": {
  "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
},
"effectiveDateTime": "2009-08-19T19:33:18-04:00",
"issued": "2009-08-19T19:33:18.715-04:00",
"dataAbsentReason": {
  "coding": [
    {
      "system": "http://terminology.hl7.org/CodeSystem/data-absent-reason",
      "code": "unknown",
      "display": "Unknown"
    }
  ],
  "text": "Unknown"
}
},
"search": {
  "mode": "match",
  "score": 1
}
},
{
  "id": "17e81f0fd34-6fc95b32-8346-4aec-9de3-072d12c992a5",
  "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Observation/17e81f0fd34-6fc95b32-8346-4aec-9de3-072d12c992a5",
  "resource": {
    "resourceType": "Observation",
    "id": "17e81f0fd34-6fc95b32-8346-4aec-9de3-072d12c992a5",

```

```

"meta": {
  "versionId": "1",
  "lastUpdated": "2022-01-22T13:19:21.396Z",
  "source": "#iLFSV7OLv0KF8dmQ",
  "profile": [
    "http://hl7.org/fhir/StructureDefinition/vitalsigns",
    "http://hl7.org/fhir/StructureDefinition/bp"
  ]
},
"status": "final",
"category": [
  {
    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/observation-category",
        "code": "vital-signs",
        "display": "vital-signs"
      }
    ]
  }
],
"code": {
  "coding": [
    {
      "system": "http://loinc.org",
      "code": "85354-9",
      "display": "Blood Pressure"
    }
  ],
  "text": "Blood Pressure"
},
"subject": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882dfb2be4b"
},

```

```

"encounter": {
  "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
},
"effectiveDateTime": "1994-05-19T19:33:18-04:00",
"issued": "1994-05-19T19:33:18.715-04:00",
"dataAbsentReason": {
  "coding": [
    {
      "system": "http://terminology.hl7.org/CodeSystem/data-absent-reason",
      "code": "unknown",
      "display": "Unknown"
    }
  ],
  "text": "Unknown"
},
"component": [
  {
    "code": {
      "coding": [
        {
          "system": "http://loinc.org",
          "code": "8462-4",
          "display": "Diastolic Blood Pressure"
        }
      ],
      "text": "Diastolic Blood Pressure"
    }
  ],
  "dataAbsentReason": {
    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/data-absent-reason",
        "code": "unknown",
        "display": "Unknown"
      }
    ]
  }
}

```

```

    ],
    "text": "Unknown"
  }
},
{
  "code": {
    "coding": [
      {
        "system": "http://loinc.org",
        "code": "8480-6",
        "display": "Systolic Blood Pressure"
      }
    ],
    "text": "Systolic Blood Pressure"
  },
  "dataAbsentReason": {
    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/data-absent-reason",
        "code": "unknown",
        "display": "Unknown"
      }
    ],
    "text": "Unknown"
  }
}
]
},
"search": {
  "mode": "match",
  "score": 1
}
},
{

```

```

    "id": "17e8502f180-3c7e4770-3630-4654-922a-3f802c2939b1",
    "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Observation/17e8502f180-3c7e4770-3630-4654-922a-3f802c2939b1",
    "resource": {
      "resourceType": "Observation",
      "id": "17e8502f180-3c7e4770-3630-4654-922a-3f802c2939b1",
      "meta": {
        "versionId": "1",
        "lastUpdated": "2022-01-23T03:37:49.696Z",
        "source": "#iLFSV7OLv0KF8dmQ",
        "profile": [
          "http://hl7.org/fhir/StructureDefinition/vitalsigns",
          "http://hl7.org/fhir/StructureDefinition/hearttrate"
        ]
      },
      "status": "final",
      "category": [
        {
          "coding": [
            {
              "system": "http://terminology.hl7.org/CodeSystem/observation-category",
              "code": "vital-signs",
              "display": "vital-signs"
            }
          ]
        }
      ],
      "code": {
        "coding": [
          {
            "system": "http://loinc.org",
            "code": "8867-4",
            "display": "Heart rate"
          }
        ]
      },

```

```

    "text": "Heart rate"
  },
  "subject": {
    "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
  },
  "encounter": {
    "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
  },
  "effectiveDateTime": "2003-07-10T19:33:18-04:00",
  "issued": "2003-07-10T19:33:18.715-04:00",
  "dataAbsentReason": {
    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/data-absent-reason",
        "code": "unknown",
        "display": "Unknown"
      }
    ],
    "text": "Unknown"
  }
},
"search": {
  "mode": "match",
  "score": 1
}
]
}

```

Request: Organization

GET [https://\[Base_URL\]/Organization/\[id\]](https://[Base_URL]/Organization/[id])

e.g. [https://\[Base_URL\]/Organization/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c](https://[Base_URL]/Organization/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c)

Response

The resource returned is expected to conform to the US Core Organization Profile

Code: 200 (Success)

```
{
  "resourceType": "Organization",
  "id": "1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c",
  "meta": {
    "versionId": "2",
    "lastUpdated": "2024-11-19T12:34:19.059+05:30",
    "source": "#OtUNHgljdSif9gwC",
    "profile": [
      "http://hl7.org/fhir/us/core/StructureDefinition/us-core-organization"
    ]
  },
  "identifier": [
    {
      "system": "https://github.com/synthetichealth/synthea",
      "value": "c7459800-e42b-3144-a934-59585c9eacff"
    },
    {
      "system": "http://hl7.org/fhir/sid/us-npi",
      "value": "1144221847"
    }
  ],
  "active": true,
  "type": [
    {
      "coding": [
        {
          "system": "http://terminology.hl7.org/CodeSystem/organization-type",
          "code": "prov",
          "display": "Healthcare Provider"
        }
      ]
    }
  ]
}
```

```

    ],
    "text": "Healthcare Provider"
  }
],
"name": "PCP170967",
"telecom": [
  {
    "system": "phone",
    "value": "413-221-0546"
  }
],
"address": [
  {
    "line": [
      "1111 ELM ST"
    ],
    "city": "WEST SPRINGFIELD",
    "state": "MA",
    "postalCode": "01089-1540",
    "country": "US"
  }
]
}

```

Request: DocumentReference

GET [https://\[Base_URL\]/DocumentReference/\[id\]](https://[Base_URL]/DocumentReference/[id])

e.g. [https://\[Base_URL\]/DocumentReference/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c](https://[Base_URL]/DocumentReference/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c)

e.g. [https://\[Base_URL\]/DocumentReference?patient=id](https://[Base_URL]/DocumentReference?patient=id)

Response

The resource returned is expected to conform to the US Core DocumentReference Profile

Code: 200 (Success)

```
{
  "resourceType": "Bundle",
  "id": "09f9423b-44a7-49c0-87f7-92737aef20c5",
  "meta": {
    "lastUpdated": "2025-11-01T17:02:26.375+05:30"
  },
  "type": "searchset",
  "total": 1,
  "link": [
    {
      "relation": "self",
      "url": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/DocumentReference?_id=1811498a297-eb9f1a18-fcac-4b00-8300-3c94af0882fe"
    }
  ],
  "entry": [
    {
      "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/DocumentReference/1811498a297-eb9f1a18-fcac-4b00-8300-3c94af0882fe",
      "resource": {
        "resourceType": "DocumentReference",
        "id": "1811498a297-eb9f1a18-fcac-4b00-8300-3c94af0882fe",
        "meta": {
          "versionId": "2",
          "lastUpdated": "2025-10-31T19:24:59.701+05:30",
          "source": "#CTLPNvIMT2G1aloz",
          "profile": [
            "http://hl7.org/fhir/us/core/StructureDefinition/us-core-documentreference"
          ]
        },
        "identifier": [
          {
            "system": "urn:ietf:rfc:3986",
```

```

    "value": "urn:uuid:18114577bc3-0a8379b4-a2f5-4cf6-b637-3301e3c0b523"
  }
],
"status": "superseded",
"type": {
  "coding": [
    {
      "system": "http://loinc.org",
      "code": "11506-3",
      "display": "Progress note"
    }
  ],
  "text": "Progress note"
},
"category": [
  {
    "coding": [
      {
        "system": "http://hl7.org/fhir/us/core/CodeSystem/us-core-documentreference-category",
        "code": "clinical-note",
        "display": "Clinical Note"
      }
    ]
  }
],
"subject": {
  "reference": "Patient/18114240088-12a390c4-2550-4130-a7d2-ff4f30bb24d0"
},
"date": "1979-02-22T18:34:18.715-05:00",
"author": [
  {
    "reference": "Practitioner/18114325cd4-68093d66-195e-4b9e-8899-44bab4231e86",
    "display": "Dr. Blossom971 Christiansen251"
  }
]

```

```

],
"custodian": {
  "reference": "Organization/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c",
  "display": "PCP170967"
},
"content": [
  {
    "attachment": {
      "contentType": "text/plain",
      "data":
"CjE5NzktMDItMjIKCiMgQ2hpZWYgQ29tcGxhaW50Ci0gRGVjcmVhc2VklGlulEp1ZGdlbWVudAotIENvbmZ1
c2lrbgotIE5hc2FsIENvbmdlc3Rpb24KLSBTbmVlemluZyBGaXRzCi0gTmFzYWwgRGZlY2hhcmdlCi0gQ291Z2g
KCgojIEhpc3Rvcnkgb2YgUHJlc2VudCBJbGxuZXNzCkR1c3RpbjMxIGlzlGEgMzggeWVhci1vbGQgbm9uLWhp
c3BhbmJlIHdoaXRlIG1hbGUuIFBhdGllbnQgaGFzIGlzlGEgaGlzdG9yeSBvZiBhY3V0ZSBhbGxlcmdpYyByZWJdGlv
biwgYWN1dGUgYnJvbmNoaXRpcyAoZGlzb3JkZXIpLCB5YWNlcmF0aW9uIG9mIHROaWdoLCBvdGI0aXMgb
WVkaWEsiH2pcmFsiHNpbnVzaXRpcyAoZGlzb3JkZXIpLCBzcHJhaW4gb2YgYW5rbGUslGZpcnN0IGRIZ3JIZSB
idXJuLCBmcmFjdHVyZSBvZiBhbmtsZSwgY2hpbGRob29klGFzdGhtYS4KCIMgU29jaWFsIEhpc3RvcnkKUGF0a
WVudCBpcyBtYXJyaWVvLiBQYXRpZW50IGlzlHRlbXBvcmlFyaWx5IGhvbWVsZXNzLiBQYXRpZW50IGlzlGFuIG
FjdGIZ2ZSBzbW9rZXIgaW5klGlzlGFuIGFsY29ob2xpYy4gUGF0aWVudCBpZGVudGlmaWVzIGFzIGhldGVyb3NI
eHVhbC4KCIBhdGllbnQgY29tZXMgZnJvbSBhIG1pZGRsZSBzb2Npb2Vjb25vbWljIGJhY2tncm91bmQulFBhdG
llbnQgZGlzIG5vdCBmaW5pc2ggaGlhbnQgY2hvb2wulFBhdGllbnQgY3VycmVudGx5IGhscyBBZXRuYS4KCIM
gQWxsZXJnaWVzCmxhdGV4IGFsbGVyZ3ksIGFsbGVyZ3kgdG8gc295YQoKlyBNZWRpY2F0aW9ucwphY2V0Y
W1pbm9waGVudMyNSBtZyBvcmlHRhYmxldDsgMTIwIGFjdHVhdCBmbHV0aWNhc29uZSBwcm9waW9u
YXRlIDAuMDQ0IG1nL2FjdHVhdCBtZXRIcmVklGRvc2UgaW5oYWxlcsjsgbmRhMDIwNTAzIDlwMCBhY3R1YXQ
gYWxidXRlcm9sIDAuMDkgbWcvYWN0dWF0IG1ldGVyZWQgZG9zZSBpbmhhbGVyOyBwcmVkbmlzb25IIDU
gbWcgb3JhbCB0YWJsZXQ7IGFzcGlyYW4gODEgbWcgb3JhbCB0YWJsZXQKCIMgQXNzZXNzbWVudCBhbm
QgUGxhbgoKCIMjIFBsYW4KCIRoZSBmb2xsb3dpbmcmcHJvY2VkdXJlcyB3ZXJlIGNvbmR1Y3RIZDoKLSBtZWR
pY2F0aW9uIHJlY29uY2lsaWF0aW9uLChwcm9jZWR1cmUpCIRoZSBwYXRpZW50IHdhcyBwcmVzY3JpYmVklH
RoZSBmb2xsb3dpbmcmcgbWVkaWNhdGlvbnM6Ci0gMTIwIGFjdHVhdCBmbHV0aWNhc29uZSBwcm9waW9uY
XRlIDAuMDQ0IG1nL2FjdHVhdCBtZXRIcmVklGRvc2UgaW5oYWxlcsjsgotIG5kYTAYMDUwMyAyMDAgYWN0d
WF0IGFsYnV0ZXJvbCAwLjA5IG1nL2FjdHVhdCBtZXRIcmVklGRvc2UgaW5oYWxlcsjsgo="
    },
    "format": {
      "system": "http://ihe.net/fhir/ValueSet/IHE.FormatCode.codesystem",
      "code": "urn:ihe:iti:xds:2017:mimeTypeSufficient",
      "display": "mimeType Sufficient"
    }
  }
],
"context": {
  "encounter": [

```

```

{
  "reference": "Encounter/1811439a512-5e32771c-2a1f-4e1d-a42d-db93aba8556d"
}
],
"period": {
  "start": "1979-02-22T18:33:18-05:00",
  "end": "1979-02-22T18:48:18-05:00"
}
}
},
"search": {
  "mode": "match"
}
}
]
}

```

Request: ServiceRequest

This is a FHIR resource that takes care of the following CCD sections:

1. Reason For Referral
2. SDOH Intervention

GET [https://\[Base_URL\]/ServiceRequest/\[id\]](https://[Base_URL]/ServiceRequest/[id])

e.g. [https://\[Base_URL\]/ServiceRequest/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c](https://[Base_URL]/ServiceRequest/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c)

e.g. [https://\[Base_URL\]/ServiceRequest?patient=id](https://[Base_URL]/ServiceRequest?patient=id)

Response

The resource returned is expected to conform to the US Core ServiceRequest Profile

Code: 200 (Success)

```

{
  "resourceType": "Bundle",
  "id": "cc11ac6e-4c04-4545-a83f-4566b9160ada",
  "meta": {

```

```

    "lastUpdated": "2025-11-01T17:09:56.724+05:30"
  },
  "type": "searchset",
  "total": 1,
  "link": [
    {
      "relation": "self",
      "url": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/ServiceRequest?_id=18f0b6042c5-d10448a0-792e-45ac-900c-fb8b9504a256"
    }
  ],
  "entry": [
    {
      "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/ServiceRequest/18f0b6042c5-d10448a0-792e-45ac-900c-fb8b9504a256",
      "resource": {
        "resourceType": "ServiceRequest",
        "id": "18f0b6042c5-d10448a0-792e-45ac-900c-fb8b9504a256",
        "meta": {
          "versionId": "1",
          "lastUpdated": "2024-11-19T12:34:23.708+05:30",
          "source": "#OtUNHgljdSif9gwC",
          "profile": [
            "http://hl7.org/fhir/us/core/StructureDefinition/us-core-servicerequest"
          ]
        },
        "identifier": [
          {
            "system": "https://github.com/inferno-framework/us-core-test-kit",
            "value": "ba89c2af-ce2b-405d-8eff-793eabe1c38d"
          }
        ],
        "status": "completed",
        "intent": "order",
        "category": [

```

```

{
  "coding": [
    {
      "system": "http://snomed.info/sct",
      "code": "387713003",
      "display": "Surgical procedure"
    }
  ],
  "text": "Surgical procedure"
},
{
  "code": {
    "coding": [
      {
        "system": "http://snomed.info/sct",
        "code": "233169004",
        "display": "Automatic defibrillator procedure (procedure)"
      }
    ]
  },
  "subject": {
    "reference": "Patient/18114240088-12a390c4-2550-4130-a7d2-ff4f30bb24d0"
  },
  "encounter": {
    "reference": "Encounter/1811439a512-5e32771c-2a1f-4e1d-a42d-db93aba8556d"
  },
  "occurrencePeriod": {
    "start": "2021-11-20",
    "end": "2021-12-04"
  },
  "authoredOn": "2015-03-30",
  "requester": {
    "reference": "Practitioner/18114325cd4-68093d66-195e-4b9e-8899-44bab4231e86"
  },

```

```

"reasonCode": [
  {
    "coding": [
      {
        "system": "http://snomed.info/sct",
        "code": "71908006",
        "display": "Ventricular fibrillation (disorder)"
      }
    ],
    "text": "Bradycardia"
  }
],
"search": {
  "mode": "match"
}
]
}

```

Request: Coverage

GET [https://\[Base_URL\]/Coverage/\[id\]](https://[Base_URL]/Coverage/[id])

e.g. [https://\[Base_URL\]/Coverage/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c](https://[Base_URL]/Coverage/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c)

e.g. [https://\[Base_URL\]/Coverage?patient=id](https://[Base_URL]/Coverage?patient=id)

Response

The resource returned is expected to conform to the US Core Coverage Profile

Code: 200 (Success)

```

{
  "resourceType": "Bundle",
  "id": "f2ecd6e3-1826-4af9-8bef-beded3c61282",

```

```

"meta": {
  "lastUpdated": "2025-11-01T17:01:39.511+05:30"
},
"type": "searchset",
"total": 1,
"link": [
  {
    "relation": "self",
    "url": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Coverage/_search"
  }
],
"entry": [
  {
    "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Coverage/18f0a59193c-e99cbeb5-f2c3-406c-9f6d-a84ce586ff9a",
    "resource": {
      "resourceType": "Coverage",
      "id": "18f0a59193c-e99cbeb5-f2c3-406c-9f6d-a84ce586ff9a",
      "meta": {
        "versionId": "1",
        "lastUpdated": "2024-11-19T12:34:21.915+05:30",
        "source": "#OtUNHgljdSif9gwC",
        "profile": [
          "http://hl7.org/fhir/us/core/StructureDefinition/us-core-coverage"
        ]
      },
      "identifier": [
        {
          "type": {

```

```

    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/v2-0203",
        "code": "MB"
      }
    ],
    "system": "https://github.com/inferno-framework/us-core-test-kit",
    "value": "f4a375d2-4e53-4f81-ba95-345e7573b550"
  }
],
"status": "active",
"type": {
  "coding": [
    {
      "system": "https://nahdo.org/sopt",
      "code": "3712",
      "display": "PPO"
    },
    {
      "system": "http://terminology.hl7.org/CodeSystem/v3-ActCode",
      "code": "PPO",
      "display": "preferred provider organization policy"
    }
  ],
  "text": "PPO"
},
"subscriberId": "888009335",
"beneficiary": {

```

```
"reference": "Patient/18114240088-12a390c4-2550-4130-a7d2-ff4f30bb24d0"
},
"relationship": {
  "coding": [
    {
      "system": "http://terminology.hl7.org/CodeSystem/subscriber-relationship",
      "code": "self"
    }
  ],
  "text": "Self"
},
"period": {
  "start": "2020-01-01"
},
"payor": [
  {
    "reference": "Organization/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c",
    "display": "PCP87052"
  }
],
"class": [
  {
    "type": {
      "coding": [
        {
          "system": "http://terminology.hl7.org/CodeSystem/coverage-class",
          "code": "group",
          "display": "Group"
        }
      ]
    }
  }
]
```

```

    ],
    "text": "An employee group"
  },
  "value": "xyz",
  "name": "XYZ employee Group Plan"
},
{
  "type": {
    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/coverage-class",
        "code": "plan",
        "display": "Plan"
      }
    ],
    "text": "A specific suite of benefits."
  },
  "value": "gold",
  "name": "Acme Gold Plan"
}
]
},
"search": {
  "mode": "match"
}
}
]
}

```

Request: MedicationDispense

GET [https://\[Base_URL\]/MedicationDispense/\[id\]](https://[Base_URL]/MedicationDispense/[id])

e.g. [https://\[Base_URL\]/MedicationDispense/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c](https://[Base_URL]/MedicationDispense/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c)

e.g. [https://\[Base_URL\]/MedicationDispense?patient=id](https://[Base_URL]/MedicationDispense?patient=id)

Response

The resource returned is expected to conform to the US Core MedicationDispense Profile

Code: 200 (Success)

```
{
  "resourceType": "Bundle",
  "id": "0ef99c93-a771-48f6-a2a2-77fb3d0ef13a",
  "meta": {
    "lastUpdated": "2025-11-01T17:02:56.829+05:30"
  },
  "type": "searchset",
  "total": 1,
  "link": [
    {
      "relation": "self",
      "url": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/MedicationDispense?patient=18114240088-12a390c4-2550-4130-a7d2-ff4f30bb24d0"
    }
  ],
  "entry": [
    {
      "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/MedicationDispense/3f4923e1-109a-4af2-93e2-dd744b6dcd4c",
      "resource": {
        "resourceType": "MedicationDispense",
```

```
"id": "3f4923e1-109a-4af2-93e2-dd744b6dcd4c",
"meta": {
  "versionId": "3",
  "lastUpdated": "2024-12-09T11:08:54.204+05:30",
  "source": "#rnoQ8mqRJg2Z0fqP",
  "profile": [
    "http://hl7.org/fhir/us/core/StructureDefinition/us-core-medicationdispense"
  ]
},
"identifier": [
  {
    "system": "https://github.com/inferno-framework/us-core-test-kit",
    "value": "1bc16d93-9ac8-4dc7-ad90-e649c01f8b10"
  }
],
"status": "completed",
"medicationReference": {
  "reference": "Medication/18114a0741e-2ab42c88-6cd0-49e2-8147-fb9954e1f4e1",
  "display": "10 ML Furosemide 10 MG/ML Injection"
},
"subject": {
  "reference": "Patient/18114240088-12a390c4-2550-4130-a7d2-ff4f30bb24d0"
},
"performer": [
  {
    "actor": {
      "reference": "Practitioner/18114325cd4-68093d66-195e-4b9e-8899-44bab4231e86",
      "display": "Dr. Willena258 Oberbrunner298"
    }
  }
]
```

```

    }
  ],
  "authorizingPrescription": [
    {
      "reference": "MedicationRequest/18114a6a311-59175c8a-f101-4638-953c-82f783333830",
      "display": "10 ML Furosemide 10 MG/ML Injection"
    }
  ],
  "type": {
    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/v3-ActCode",
        "code": "FFC",
        "display": "First Fill - Complete"
      }
    ]
  },
  "quantity": {
    "value": 480,
    "unit": "mL",
    "system": "http://unitsofmeasure.org",
    "code": "mL"
  },
  "whenHandedOver": "2016-12-08T06:38:52Z",
  "dosageInstruction": [
    {
      "text": "10 mL bid",
      "timing": {
        "repeat": {

```

```

    "boundsPeriod": {
      "start": "2008-04-05"
    },
    "frequency": 2,
    "period": 1,
    "periodUnit": "d"
  }
},
"doseAndRate": [
  {
    "doseQuantity": {
      "value": 10,
      "unit": "ml",
      "system": "http://unitsofmeasure.org",
      "code": "mL"
    }
  }
]
}
]
},
"search": {
  "mode": "match"
}
}
]
}

```

Request: Specimen

GET [https://\[Base_URL\]/Specimen/\[id\]](https://[Base_URL]/Specimen/[id])

e.g. [https://\[Base_URL\]/Specimen/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c](https://[Base_URL]/Specimen/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c)

Response

The resource returned is expected to conform to the US Core SpecimenProfile

Code: 200 (Success)

```
{
  "resourceType": "Specimen",
  "id": "18f0abe1999-9afa3deb-e111-4a8f-a8c3-914d100b49d8",
  "meta": {
    "versionId": "1",
    "lastUpdated": "2024-11-19T12:34:22.045+05:30",
    "source": "#OtUNHgljdSif9gwC",
    "profile": [
      "http://hl7.org/fhir/us/core/StructureDefinition/us-core-specimen"
    ]
  },
  "identifier": [
    {
      "system": "https://github.com/inferno-framework/us-core-test-kit",
      "value": "82115a91-04c9-4654-acb5-72f2806e00fc"
    }
  ],
  "accessionIdentifier": {
    "system": "https://github.com/inferno-framework/us-core-test-kit",
    "value": "7446d969-3293-4f51-8d48-034121156b4e"
  },
  "type": {
    "coding": [
      {
        "system": "http://snomed.info/sct",
        "code": "122575003",
```

```

      "display": "Urine specimen (specimen)"
    }
  ]
},
"subject": {
  "reference": "Patient/18114240088-12a390c4-2550-4130-a7d2-ff4f30bb24d0"
},
"collection": {
  "collectedDateTime": "2005-07-05T18:07:12Z",
  "bodySite": {
    "coding": [
      {
        "system": "http://snomed.info/sct",
        "code": "106004",
        "display": "Posterior carpal region"
      }
    ]
  }
},
"container": [
  {
    "type": {
      "coding": [
        {
          "system": "http://snomed.info/sct",
          "code": "702244006",
          "display": "Sterile urine specimen container"
        }
      ]
    },
    "text": "Sterile urine specimen container"
  }
],
"condition": [

```

```
{
  "coding": [
    {
      "system": "http://terminology.hl7.org/CodeSystem/v2-0493",
      "code": "AUT",
      "display": "Autolyzed"
    }
  ]
}
```

Request: RelatedPerson

GET [https://\[Base_URL\]/RelatedPerson/\[id\]](https://[Base_URL]/RelatedPerson/[id])

e.g. [https://\[Base_URL\]/RelatedPerson/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c](https://[Base_URL]/RelatedPerson/1811435a7ea-7a1b883e-e0f0-4cb2-b938-4b6fd8d60f0c)

Response

The resource returned is expected to conform to the US Core RelatedPerson Profile

Code: 200 (Success)

```
{
  "resourceType": "RelatedPerson",
  "id": "re1a0ab5-5e30-46b2-868d-7e56c252b1f0",
  "meta": {
    "versionId": "1",
    "lastUpdated": "2024-12-09T12:40:19.084+05:30",
    "source": "#sCPUceOmdiu7NGzZ",
    "profile": [
      "http://hl7.org/fhir/us/core/StructureDefinition/us-core-relatedperson"
    ]
  },
  "active": false,
  "patient": {
```

```
"reference": "Patient/18114240088-12a390c4-2550-4130-a7d2-ff4f30bb24d0"
},
"relationship": [
  {
    "coding": [
      {
        "system": "http://terminology.hl7.org/CodeSystem/v3-RoleCode",
        "code": "ROOM",
        "display": "Roommate"
      }
    ]
  }
],
"name": [
  {
    "use": "official",
    "family": "Jefferson174",
    "given": [
      "Ronald408",
      "MacGyver246"
    ],
    "prefix": [
      "Mr."
    ]
  }
],
"telecom": [
  {
    "system": "phone",
    "value": "555-164-3435",
    "use": "home"
  }
],
"address": [
```

```

{
  "extension": [
    {
      "url": "http://hl7.org/fhir/StructureDefinition/geolocation",
      "extension": [
        {
          "url": "latitude",
          "valueDecimal": 42.33734137767395
        },
        {
          "url": "longitude",
          "valueDecimal": -71.0633497363917
        }
      ]
    }
  ],
  "line": [
    "445 Braun Mews"
  ],
  "city": "Boston",
  "state": "MA",
  "postalCode": "02110",
  "country": "US",
  "period": {
    "start": "2019-02-07"
  }
}

```

Request: Observation

This is a FHIR resource that takes care of the following CCDA sections:

1. Pregnancy Status
2. Pregnancy Intent
3. Occupation

4. Respiratory Rate
5. Heart Rate
6. Body Temperature
7. Pediatric Weight For Height
8. Pulse Oximetry
9. Head Circumference
10. Body Height
11. BMI
12. Screening Assessment
13. Blood Pressure
14. Clinical Result
15. Pediatric BMI
16. Pediatric Head Occipital Frontal Circumference Percentile
17. Observation Body Weight
18. SDOH Assessment
19. Basic Occupation and Occupation Industry
20. Tribal Affiliation
21. Functional Status
22. Cognitive Status
23. Disability Status
24. Gender Identity
25. Sexual Orientation
26. Sex at Birth

In the below GET URL, you need to mention appropriate [code] and you will get the required bundle.

GET [https://\[Base_URL\]/Observation?code=\[code\]&patient=\[id\]](https://[Base_URL]/Observation?code=[code]&patient=[id])

e.g. [https://\[Base_URL\]/Observation?code=72166-2&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b](https://[Base_URL]/Observation?code=72166-2&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b)

Response

Each resource returned from the first search is expected to conform to the US Core [Smoking Status Observation](#) Profile

Code: 200 (Success)

```
{
  "resourceType": "Bundle",
  "id": "2ea85043-b450-40f1-8d2d-751af1c10abb",
  "type": "searchset",
  "total": 2,
  "link": [
    {
```

```

    "relation": "self",
    "url": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Observation?_count=10&code=72166-2&patient=17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b&_page=1"
  }
],
"entry": [
  {
    "id": "17e77809e90-d5cfd26d-ab5d-4f36-bc32-1ee799b2c278",
    "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Observation/17e77809e90-d5cfd26d-ab5d-4f36-bc32-1ee799b2c278",
    "resource": {
      "resourceType": "Observation",
      "id": "17e77809e90-d5cfd26d-ab5d-4f36-bc32-1ee799b2c278",
      "meta": {
        "versionId": "1",
        "lastUpdated": "2022-01-20T12:40:24.976Z",
        "profile": [
          "http://hl7.org/fhir/us/core/StructureDefinition/us-core-smokingstatus"
        ]
      },
      "status": "final",
      "category": [
        {
          "coding": [
            {
              "system": "http://terminology.hl7.org/CodeSystem/observation-category",
              "code": "survey",
              "display": "survey"
            }
          ]
        }
      ],
      "code": {
        "coding": [
          {

```

```

    "system": "http://loinc.org",
    "code": "72166-2",
    "display": "Tobacco smoking status NHIS"
  }
],
"text": "Tobacco smoking status NHIS"
},
"subject": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
},
"encounter": {
  "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
},
"effectivePeriod": {
  "start": "1940-09-05T19:33:18-04:00",
  "end": "1940-09-05T20:33:18-04:00"
},
"issued": "1940-09-05T19:33:18.715-04:00",
"valueCodeableConcept": {
  "coding": [
    {
      "system": "http://snomed.info/sct",
      "code": "266919005",
      "display": "Never smoker"
    }
  ],
  "text": "Never smoker"
}
},
"search": {
  "mode": "match",
  "score": 1
}
},

```

```

{
  "id": "17e7785ea69-6c9f88d9-0553-4216-be17-f908ebe0b371",
  "fullUrl": "https://fhir.iPatientCare.net:9443/fhirserver/fhir/Observation/17e7785ea69-6c9f88d9-0553-4216-be17-f908ebe0b371",
  "resource": {
    "resourceType": "Observation",
    "id": "17e7785ea69-6c9f88d9-0553-4216-be17-f908ebe0b371",
    "meta": {
      "versionId": "1",
      "lastUpdated": "2022-01-20T12:46:12.073Z",
      "profile": [
        "http://hl7.org/fhir/us/core/StructureDefinition/us-core-smokingstatus"
      ]
    },
    "status": "final",
    "category": [
      {
        "coding": [
          {
            "system": "http://terminology.hl7.org/CodeSystem/observation-category",
            "code": "survey",
            "display": "survey"
          }
        ]
      }
    ],
    "code": {
      "coding": [
        {
          "system": "http://loinc.org",
          "code": "72166-2",
          "display": "Tobacco smoking status NHIS"
        }
      ]
    },
    "text": "Tobacco smoking status NHIS"
  }
}

```

```

},
"subject": {
  "reference": "Patient/17e77411359-8837d878-ea1e-45b5-bb89-882fdb2be4b"
},
"encounter": {
  "reference": "Encounter/17e77566282-0c9c236b-f0bd-45e0-9c61-80991a376967"
},
"effectivePeriod": {
  "start": "1940-09-05T19:33:18-04:00",
  "end": "1940-09-05T20:33:18-04:00"
},
"issued": "1940-09-05T19:33:18.715-04:00",
"valueCodeableConcept": {
  "coding": [
    {
      "system": "http://snomed.info/sct",
      "code": "266919005",
      "display": "Never smoker"
    }
  ],
  "text": "Never smoker"
}
},
"search": {
  "mode": "match",
  "score": 1
}
}
]
}

```

Profile audience and scope

This profile is intended to be used by developers of apps that need to access user identity information or other FHIR resources by requesting authorization from OAuth 2.0 compliant authorization servers.

The profile defines a method through which an app requests authorization to access a FHIR resource, and then uses that authorization to retrieve the resource. Synchronization of patient context is not addressed; for use cases that require context synchronization (e.g., learning about when the in-context patient changes within an EHR session) In other words, if the patient chart is changed during the session, the application will not inherently be updated.

Security and Privacy Considerations

App Protection

The app is responsible for protecting itself from potential misbehaving or malicious values passed to its redirect URL (e.g., values injected with executable code, such as SQL) and for protecting authorization codes, access tokens, and refresh tokens from unauthorized access and use. The app developer must be aware of potential threats, such as malicious apps running on the same platform, counterfeit authorization servers, and counterfeit resource servers, and implement countermeasures to help protect both the app itself and any sensitive information it may hold. For background, see the [OAuth 2.0 Threat Model and Security Considerations](#).

Specific requirements are:

- Apps SHALL ensure that sensitive information (authentication secrets, authorization codes, tokens) is transmitted ONLY to authenticated servers, over TLS-secured channels. We support TLS 1.2. Ensure that TLS 1.2 is enabled on our web server's operating system by modifying registry settings
- Apps SHALL generate an unpredictable **state** parameter for each user session; SHALL include **state** with all authorization requests; and SHALL validate the **state** value for any request sent to its redirect URL.
- An app SHALL NOT execute untrusted user-supplied inputs as code.
- App SHALL NOT forward values passed back to its redirect URL to any other arbitrary or user-provided URL (a practice known as an "open redirector").
- An app SHALL NOT store bearer tokens in cookies that are transmitted as clear text.
- Apps SHOULD persist tokens and other sensitive data in app-specific storage locations only, and SHOULD NOT persist them in system-wide-discoverable locations.
- Secure issuance of refresh token: Native applications call our OAuth api to request tokens. Since all our endpoints are on https they are secure.
- We support PKCE

Terms of Use

1. TERMS ACCEPTANCE AND REPRESENTATION

- a. Accepting the Terms. These Terms of Use ("Terms") govern your access to and use of the Certified iPatientCare EHR API, documentation, services, etc. By accessing or using the Certified iPatientCare APIs, you agree to be bound by these Terms. "Certified iPatientCare EHR API" means the API provided by iPatientCare to allow authorized access to query our Client(s) Electronic Health Record system. You represent and warrant that you are at least 18 years of age and that you possess the legal right and ability to agree to these Terms and to use the Certified iPatientCare APIs in accordance with these Terms.
- b. Entity Representation. If you are using the Certified iPatientCare API on behalf of a legal entity (i.e. a Clinical software services company), you represent that you have proper authority to act on behalf of and bind the entity to these Terms, and by accepting, you accept on behalf of the entity (and all references to "you" in the Terms refer to the entity).

2. REPRESENTATIONS AND RESPONSIBILITIES

- a. Compliance. You agree to be financially responsible for your use of the Certified iPatientCare EHR APIs and to comply with your responsibilities and obligations as stated in these Terms. You agree to comply at all times with all applicable laws, rules and regulations relating to the use of the Certified iPatientCare EHR APIs. You hereby grant iPatientCare the right to monitor and periodically audit in a reasonable manner your use of the Certified iPatientCare APIs, your App and other activities related to your obligations under these Terms.
- b. Virus Warranty. You warrant that your Apps will not contain any viruses or other malicious computer instructions, devices, or techniques that can or were designed to threaten, infect, damage, disable, or shut down the iPatientCare EHR APIs, any technology, software, solution, equipment or any computer system.

3. GENERAL

- a. Changes. iPatientCare may, in its sole and absolute discretion, make changes, modifications or updates to the Certified iPatientCare EHR API (including without limitation changes to the capabilities and tech specs), without notice to you.
- b. Global Availability. iPatientCare makes no representations that the Certified iPatientCare EHR APIs are appropriate or available for use in locations outside of the United States, and access to them from such territories is at your own risk. Those who choose to access the Certified iPatientCare APIs from locations outside of the United States do so at their own initiative and are responsible for compliance with applicable local laws.
- c. Intellectual Property Rights. You acknowledge and agree that the Certified iPatientCare EHR APIs and iPatientCare's software, products and services are proprietary in nature, that iPatientCare claims all intellectual property rights therein as well as in all modifications, enhancements and alterations thereto, and that iPatientCare neither grants nor otherwise transfers any rights of ownership therein to you or any third party. No rights or licenses are granted by iPatientCare other than those rights expressly granted in these Terms, and iPatientCare reserves all rights not expressly granted.
- d. Waiver, Release and Limitation of Liability. You acknowledge and agree that neither iPatientCare nor any of its affiliates will have any responsibility or liability with respect to your or any third party's distribution, implementation, commercialization, use, or other

form of exploitation of any app, software, solution, service or other technology created by you or another third party. As between you and iPatientCare, you are solely responsible and liable for all representations, warranties, support and other obligations made by you to any third party related to any app, software, solution, service or other technology created by you, including claims arising from product liability, breach of warranty, and intellectual property infringement. YOU HEREBY RELEASE AND FOREVER WAIVE ANY AND ALL CLAIMS YOU MAY HAVE AGAINST IPATIENTCARE, ITS OFFICERS, DIRECTORS, EMPLOYEES, AGENTS, INFORMATION PROVIDERS OR SUPPLIERS FOR LOSSES OR DAMAGES YOU SUSTAIN IN CONNECTION WITH YOUR USE OF THE CERTIFIED IPATIENTCARE EHR API AND ANY IPATIENTCARE WEBSITES. IPATIENTCARE MAKES NO REPRESENTATIONS ABOUT THE SUITABILITY, RELIABILITY, AVAILABILITY, TIMELINESS AND ACCURACY OF THE CERTIFIED IPATIENTCARE APIs OR OTHER INFORMATION, TECHNOLOGY, SOFTWARE, PRODUCTS AND SERVICES PROVIDED BY IPATIENTCARE FOR ANY PURPOSE. ALL SUCH TECHNOLOGY, INFORMATION, SOFTWARE, PRODUCTS AND SERVICES ARE PROVIDED "AS IS" WITHOUT WARRANTY OF ANY KIND. IPATIENTCARE HEREBY DISCLAIMS ALL WARRANTIES AND CONDITIONS WITH REGARD TO THE CERTIFIED IPATIENTCARE APIs, OTHER TECHNOLOGY, INFORMATION, SOFTWARE, SOLUTIONS, PRODUCTS AND SERVICES, INCLUDING ALL IMPLIED WARRANTIES AND CONDITIONS OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, TITLE AND NON-INFRINGEMENT. TO THE FULLEST EXTENT PERMITTED BY LAW, IN NO EVENT SHALL IPATIENTCARE, ITS OFFICERS, DIRECTORS, EMPLOYEES, AGENTS, AND/OR ITS SUPPLIERS BE LIABLE FOR ANY DIRECT, INDIRECT, PUNITIVE, INCIDENTAL, SPECIAL, CONSEQUENTIAL DAMAGES OR ANY OTHER DAMAGES WHATSOEVER INCLUDING, WITHOUT LIMITATION, DAMAGES FOR LOSS OF USE, DATA OR PROFITS, ARISING OUT OF OR IN ANY WAY CONNECTED WITH THE ACCESS, USE OR PERFORMANCE OF THE CERTIFIED IPATIENTCARE APIs, WITH THE DELAY OR INABILITY TO USE THE CERTIFIED IPATIENTCARE EHR APIs OR RELATED TECHNOLOGY, SOFTWARE OR SERVICES, THE PROVISION OF OR FAILURE TO PROVIDE THE CERTIFIED IPATIENTCARE APIs, SOFTWARE OR SERVICES, OR FOR ANY INFORMATION, SOFTWARE, PRODUCTS AND SERVICES OBTAINED FROM IPATIENTCARE, OR OTHERWISE ARISING OUT OF THE USE OF THE CERTIFIED IPATIENTCARE APIs, WHETHER BASED ON CONTRACT, TORT, NEGLIGENCE, STRICT LIABILITY OR OTHERWISE, EVEN IF IPATIENTCARE HAS BEEN ADVISED OF THE POSSIBILITY OF DAMAGES. BECAUSE SOME STATES/JURISDICTIONS DO NOT ALLOW THE EXCLUSION OR LIMITATION OF LIABILITY FOR CONSEQUENTIAL OR INCIDENTAL DAMAGES, THE ABOVE LIMITATION MAY NOT APPLY TO YOU. IF YOU ARE DISSATISFIED WITH ANY PORTION OF THE CERTIFIED CLINCMAX EHR APIs, OR WITH THESE TERMS OF USE, YOUR SOLE AND EXCLUSIVE REMEDY IS TO DISCONTINUE USING THE CERTIFIED IPATIENTCARE EHR APIs. NOTWITHSTANDING THE FOREGOING PARAGRAPH, THE TOTAL LIABILITY OF IPATIENTCARE, ITS OFFICERS, DIRECTORS, EMPLOYEES, AGENTS, INFORMATION PROVIDERS AND SUPPLIERS, IF ANY, FOR LOSSES OR DAMAGES SHALL NOT EXCEED THE FEES PAID BY YOU FOR THE USE OF THE PARTICULAR TECHNOLOGY, SOFTWARE, PRODUCT, INFORMATION OR SERVICE PROVIDED BY IPATIENTCARE.

- e. Indemnification. YOU AGREE TO INDEMNIFY, DEFEND AND HOLD HARMLESS IPATIENTCARE, ITS OFFICERS, DIRECTORS, EMPLOYEES, AGENTS, INFORMATION PROVIDERS AND SUPPLIERS FROM AND AGAINST ALL CLAIMS, LIABILITIES, LOSSES,

EXPENSES, DAMAGES AND COSTS, INCLUDING REASONABLE ATTORNEYS' FEES, RESULTING FROM ANY VIOLATION OF THESE TERMS OR ANY ACTIVITY RELATED TO YOUR USE OF THE CERTIFIED IPATIENTCARE APIS OR THE IPATIENTCARE WEB SITES.

- f. Term and Termination. Either you or iPatientCare may terminate your right to use the Certified iPatientCare EHR APIs at any time, with or without cause, upon notice. iPatientCare also reserves the right to disable your API access in a production environment at any time, with or without cause. iPatientCare reserves the right to disable access to the iPatientCare EHR APIs if your App poses any security, privacy, or patient safety risks. The provisions concerning Indemnification, Waiver, Release and Limitation of Liability, and General shall survive any termination of these Terms.
 - g. Governing Law. These Terms are governed by US federal law or the laws of the State of Ohio.
4. Application developer affirmations to Certified API Developers regarding the ability of their applications to secure a refresh token, a client secret, or both, must be treated in a good faith manner